

# UC Merced

## UC Merced Previously Published Works

### Title

In Vitro Culturing and Screening of Candida albicans Biofilms

### Permalink

<https://escholarship.org/uc/item/12776923>

### Journal

Current Protocols in Microbiology, 50(1)

### ISSN

1934-8525

### Authors

Gulati, Megha  
Lohse, Matthew B  
Ennis, Craig L  
[et al.](#)

### Publication Date

2018-08-01

### DOI

10.1002/cpmc.60

Peer reviewed

# *In Vitro* Culturing and Screening of *Candida albicans* Biofilms

Megha Gulati,<sup>1</sup> Matthew B. Lohse,<sup>2,3</sup> Craig L. Ennis,<sup>1,4</sup> Ruth E. Gonzalez,<sup>1</sup> Austin M. Perry,<sup>1,4</sup> Priyanka Bapat,<sup>1,4</sup> Ashley Valle Arevalo,<sup>1,4</sup> Diana L. Rodriguez,<sup>1,4</sup> and Clarissa J. Nobile<sup>1,5</sup>

<sup>1</sup>Department of Molecular and Cell Biology, School of Natural Sciences, University of California, Merced, Merced, California

<sup>2</sup>Department of Microbiology and Immunology, University of California, San Francisco, San Francisco, California

<sup>3</sup>Department of Biology, BioSynthesis, Inc., San Francisco, California

<sup>4</sup>Quantitative and Systems Biology Graduate Program, University of California, Merced, Merced, California

<sup>5</sup>Corresponding authors: [cnobile@ucmerced.edu](mailto:cnobile@ucmerced.edu)

*Candida albicans* is a normal member of the human microbiota that asymptotically colonizes healthy individuals, however it is also an opportunistic pathogen that can cause severe infections, especially in immunocompromised individuals. The medical impact of *C. albicans* depends, in part, on its ability to form biofilms, communities of adhered cells encased in an extracellular matrix. Biofilms can form on both biotic and abiotic surfaces, such as tissues and implanted medical devices. Once formed, biofilms are highly resistant to antifungal agents and the host immune system, and can act as a protected reservoir to seed disseminated infections. Here, we present several *in vitro* biofilm protocols, including protocols that are optimized for high-throughput screening of mutant libraries and antifungal compounds. We also present protocols to examine specific stages of biofilm development and protocols to evaluate interspecies biofilms that *C. albicans* forms with interacting microbial partners. © 2018 by John Wiley & Sons, Inc.

Keywords: biofilm protocols • biofilm methods • biofilm screens • *Candida albicans* • interspecies biofilms

## How to cite this article:

Gulati, M., Lohse, M. B., Ennis, C. L., Gonzalez, R. E., Perry, A. M., Bapat, P., Valle Arevalo, A., Rodriguez, D. L., & Nobile, C. J. (2018). *In vitro* culturing and screening of *Candida albicans* biofilms. *Current Protocols in Microbiology*, e60. doi: 10.1002/cpmc.60

## INTRODUCTION


*Candida albicans* is a normal member of the human microbiota that asymptotically colonizes several niches of the body (e.g., skin, ears, nasal cavity, mucosal membranes, gastrointestinal and urogenital tracts; Douglas, 2003; Gulati & Nobile, 2016; Nobile & Johnson, 2015). *C. albicans* is also one of the few fungal species that can cause disease in humans, which can range from superficial mucosal and dermal infections to severe disseminated bloodstream and deep-seated tissue infections (Douglas, 2003; Kim & Sudbery, 2011; Kullberg & Oude Lashof, 2002; Nobile & Johnson, 2015). These infections are especially serious in immunocompromised individuals (Calderone & Fonzi, 2001; Douglas, 2003; Lopez-Ribot, 2005). One important virulence trait of

*C. albicans* is its ability to form biofilms, structured communities of cells that are encased in an extracellular matrix and are adhered to a surface (Chandra et al., 2001; Douglas, 2002; Kumamoto, 2002; Lohse, Gulati, Johnson, & Nobile, 2018; Nobile & Johnson, 2015). *C. albicans* biofilms can form on both biotic and abiotic surfaces, such as tissues and implanted medical devices, are highly resistant to physical and chemical perturbations, and serve as protected reservoirs that can seed new biofilm infections as well as disseminated (non-biofilm) infections (Douglas, 2002, 2003; Gulati & Nobile, 2016).

*C. albicans* produces structured biofilms consisting of multiple cell types (spherical yeast-form cells, oval pseudohyphal cells, and cylindrical hyphal cells; Douglas, 2003; Gulati & Nobile, 2016). *C. albicans* biofilm formation proceeds through four distinct stages: (1) adherence, where yeast-form cells attach to a surface to seed a biofilm; (2) initiation, where the adhered cells proliferate on the surface to form an anchoring basal layer; (3) maturation, where cells filament and continue to proliferate, leading to a several hundred micron thick biofilm with layers of intercalating hyphae, pseudohyphae, and yeast-form cells encased in an extracellular matrix; and (4) dispersion, where yeast-form cells are released from the biofilm to seed new sites (Baillie & Douglas, 1999; Chandra et al., 2001; Gulati & Nobile, 2016; Lohse et al., 2018; Nobile & Johnson, 2015; Uppuluri et al., 2010).

To study biofilm formation in the lab, typical *in vitro* biofilm assays involve an adherence step where cells first adhere to a solid surface, a wash step to remove non- and weakly adhered cells, and a maturation step where the adhered cells develop into the biofilm. The final step of the assay entails some sort of measurement of the resulting biofilm (e.g., optical density measurements using a plate reader or microscopic measurements using a confocal scanning laser microscope). For the majority of *in vitro C. albicans* biofilm assays, the biofilm is exposed to either shaking conditions (using a shaking incubator) or to continuous flow across the biofilm surface (using a microfluidic device) throughout the adherence and maturation steps (Lohse et al., 2017; Tournu & Van Dijck, 2012). Specific *in vitro* biofilm assays vary in terms of how the growth of the biofilm is evaluated, such as by dry weight (Hawser & Douglas, 1994; Nobile, Andes, et al., 2006; Nobile et al., 2012; Nobile, Nett, et al., 2006), optical density (Fox et al., 2015; Lohse et al., 2017; Nobile et al., 2014; Uppuluri et al., 2010; Winter et al., 2016), cell viability (Nailis et al., 2010), or direct observations like confocal scanning laser microscopy (Nobile et al., 2012; Nobile et al., 2005). Other variables include the material the biofilm grows on, such as treated or untreated polystyrene plates (Krom & Willems, 2016; Lohse et al., 2017) or silicone squares (Nobile, Andes, et al., 2006; Nobile et al., 2012; Nobile, Nett, et al., 2006), and the specific stage of biofilm development being observed, such as adherence (Finkel et al., 2012; Lohse et al., 2017; Winter et al., 2016) or dispersion (Lohse et al., 2017; Nobile et al., 2014; Uppuluri et al., 2010). These *in vitro* biofilm assays can also be used to assess the biofilms formed by different strains, specific mutants of interest (Finkel et al., 2012; Fox et al., 2015; Nobile, Andes, et al., 2006; Nobile et al., 2012; Nobile & Mitchell, 2005; Nobile, Nett, et al., 2006; Nobile et al., 2009; Norice, Smith, Solis, Filler, & Mitchell, 2007; Richard, Nobile, Bruno, & Mitchell, 2005), or upon exposure to antifungal compounds (Lafleur et al., 2013; LaFleur, Lucumi, Napper, Diamond, & Lewis, 2011; Pierce et al., 2015; Pierce, Saville, & Lopez-Ribot, 2014).

Here we present several validated and commonly used *in vitro* biofilm protocols designed to investigate different aspects of *C. albicans* biofilm formation, each with their individual trade-offs in terms of information generated, throughput, and infrastructure requirements (Figure 1 and Table 1). On the high throughput end of the spectrum, we present an optical density-based biofilm formation assay using 96- or 384-well microtiter plates that allows for rapid high-throughput screening of large deletion libraries and



**Figure 1** Overview of the steps of a typical *in vitro* biofilm assay and its relationship to the protocols presented in this article. A typical *in vitro* biofilm assay involves an adherence step, where cells first adhere to a solid surface, a wash step to remove non- and weakly adhered cells, and a maturation step where the adhered cells develop into the mature biofilm.

testing of putative antifungal compounds. We present several variations of this assay, each designed to investigate different aspects of biofilm formation (Figure 1). We also present protocols that allow for the enumeration of live/dead cells within a biofilm, the measurement of biofilm biomass, and the collection of biofilm conditioned medium for use in downstream proteomic assays. Additional protocols presented include the formation of sexual biofilms, which are biofilms made by mating-competent cells in response to mating pheromones, and the formation of mixed-species biofilms. Finally, we also present protocols for visualizing *C. albicans* biofilms using a confocal microscope that allows for assessment of biofilm architecture, and a customizable microfluidic assay that allows for real-time visualization of biofilm formation over time in host-mimicking conditions. The collection of *in vitro* protocols outlined below are protocols that we have vetted and recommend for the assessment of specific aspects of *C. albicans* biofilm formation.

## CULTURING

In most cases, *C. albicans* strains are grown in yeast extract peptone dextrose (YEPD; also abbreviated YPD) medium at 30°C prior to biofilm assays as this condition has worked for a wide range of *C. albicans* clinical isolates and deletion strains as well as for other *Candida* species (e.g., *Candida tropicalis* and *Candida parapsilosis*). *Candida* strains of interest (wild-type and mutant strains) can be obtained from sev-

**Table 1** Protocols Presented in This Unit

Assay	Technique number	Throughput	Specialized equipment required <sup>a</sup>	Result format
Culturing	BP 1	NA		NA
Standard optical density	BP 2	High	Plate reader	OD measurement
XTT reduction	AP 1	High	Plate reader	OD measurement
Crystal violet	AP 2	High	Plate reader	OD measurement
Inhibition optical density	AP 3	High	Plate reader	OD measurement
Disruption optical density	AP 4	High	Plate reader	OD measurement
Dispersal optical density	BP 3	High	Plate reader	OD measurement
Co-culturing and analyzing with standard optical density	BP 4	High	Plate reader	OD measurement
Colony-forming units	AP 5	Medium		Colonies on plate
Silicone square	BP 5	Medium		Depends on usage
Cell adhesion	BP 6	Medium		Colonies on plate
Dry weight	BP 7	Low	Analytical balance	Weight of biofilm
Cell enumeration	BP 8	Medium		Colonies on plate
Cell enumeration using fluorescence microscopy	AP 6	Medium	Fluorescent microscope	Image
Pheromone-stimulated	BP 9	High	Plate reader	OD measurement
Microfluidics	BP 10	Medium	BioFlux EZ1000	Image/movie
Confocal imaging	BP 11	Low	Confocal microscope	Image
Confocal imaging of mixed-species biofilms	AP 7	Low	Confocal microscope	Image
Conditioned medium harvesting	BP 12	Low		Depends on usage

<sup>a</sup>Microplate shaking incubators are required for all *C. albicans* biofilm assays.

eral sources, such as the American Type Culture Collection (<https://www.atcc.org/>), the Centers for Disease Control and Prevention (CDC) Antibiotic Resistance (AR) Isolate Bank (<https://wwwn.cdc.gov/arisolatebank/Overview>), the Fungal Genetics Stock Center (<http://www.fgsc.net/>), or from various academic labs. Alternative media may be substituted for YEPD as needed, although various aspects of the protocol (e.g., duration of growth) may need to be re-optimized as a result. Likewise, strains with severe growth defects may require longer growth periods. This protocol has been optimized for a modest number of strains (i.e., one to dozens at a time); higher throughput options (e.g., overnight cultures grown in deep-well 96-well plates) may be preferable if working with larger numbers of strains at the same time.

### Materials

*Candida albicans* strain(s) of interest stored in glycerol  
 Yeast extract peptone dextrose (YEPD) liquid medium (2% Bacto™ peptone, 2% dextrose, 1% yeast extract; store at room temperature)  
 Yeast extract peptone dextrose (YEPD) agar plates (2% Bacto™ peptone, 2% dextrose, 1% yeast extract, 2% agar; store at room temperature)  
 20-ml test tubes and/or 250-ml flasks  
 Roller drum (for test tubes) and/or shaker (for flasks) in a 30°C incubator  
 Cuvettes and spectrophotometer or transparent 96-well plate and OD<sub>600</sub>-capable plate reader

1. Streak *C. albicans* strains onto YEPD plates from glycerol stocks 2 to 5 days in advance of assay. Incubate plates at 30°C, 48 hr.

*Do not use colonies that are >7 days old or store plates at 4°C as C. albicans can acquire aneuploidies under these conditions. Likewise, new plates should be streaked from glycerol stocks rather than by re-streaking cells from an existing plate.*

- 2a. For small numbers of strains/large volumes: Inoculate a single colony in 25 ml YEPD medium in a 250-ml flask from a 2- to 5-day-old plate. Grow at 30°C with shaking overnight.
- 2b. For large numbers of strains/small volume: Inoculate a single colony in 4 ml YEPD medium in a 20-ml test tube from a 2- to 5-day-old plate. Grow at 30°C with shaking or on a roller drum overnight.
3. After 16 to 18 hr, dilute overnight culture 1:20 or 1:40 and measure OD<sub>600</sub> using cuvettes and a spectrophotometer or 96-well plates and a plate reader. Dilute strains as indicated by the specific protocol being used.

*If assays will be set up >18 hr after overnight cultures were started, remove overnight cultures from the shaker or roller drum in the morning and allow to sit at room temperature until time of use (cultures should not sit for >4 hr at room temperature). Vortex cultures immediately before use. However, if possible, it is better to time the experiment such that the strains will be used prior to 18 hr of overnight growth. Depending on the assay, it may not be necessary to determine the density of the culture(s) prior to use.*

## STANDARD OPTICAL DENSITY ASSAY

The standard optical density assay offers a quick, relatively high-throughput way to screen large numbers of mutant strains for defects in biofilm formation with minimal equipment requirements (Fox et al., 2015; Lohse et al., 2017). The basic version of this assay quantifies biofilm formation based solely on the optical density of the biofilm rather than measurement of cell viability, e.g., metabolic reduction of the tetrazolium salt reagent 2,3-bis-(2-methoxy-4-nitro-5-sulphophenyl)-2H-tetrazolium-5-carboxanilide (XTT) or the uptake of dyes, such as crystal violet. This assay avoids the extra washing and dye addition steps needed for the XTT reduction or crystal violet assays while avoiding issues with physical disruption of the biofilms, incomplete reagent penetrance, or confounding issues related to the presence of metabolically inactive cells.

Variants of this assay can be used to determine the ability of compounds to prevent the formation of a biofilm or to disrupt an established biofilm. We typically use RPMI-1640 medium for this assay, however other media (e.g., Spider medium, which results in thinner biofilms, and thus can reveal subtler biofilm defects) can also be used. Although this protocol can be modified to grow biofilms on 6- and 12-well plates, using 4 ml and 2 ml of medium per well, respectively, we recommend using 96- or 384-well formats for high-throughput assays. We recommend the use of six wells per strain or condition in 96-well plate format and eight wells per strain or condition in 384-well plate format, allowing for the evaluation of 16 conditions per 96-well plate and 48 conditions per 384-well plate (including blanks and controls). This protocol has been optimized for use with *C. albicans*, however it can be adapted for use with other microbial species. When optimizing this protocol for use with other species, changes to the medium, temperature, and inoculation amounts are common conditions to vary. Deep well 96- or 384-well plates may be helpful for pre-aliquoting compound/medium mixtures, especially when performing the alternative versions of this assay.

## Materials

*Candida albicans* cultures (see Dulbecco's phosphate-buffered saline [PBS; calcium and magnesium salt free]), sterile filtered  
RPMI-1640 medium with L-glutamine and 3-(*N*-morpholino)propanesulfonic acid (MOPS), without sodium bicarbonate, pH 7.0 (see recipe)  
Transparent, sterile, flat-bottomed, non-tissue culture treated 96-well (BD Falcon, model number 351172) or 384-well (Thermo Fisher Scientific, Thermo Scientific™ 242765) microtiter plates  
Breathe-Easy® sealing membranes (Diversified Biotech BEM-1)  
Sterile reagent reservoirs (USA Scientific, cat. no. 2321-2530)  
Multi-channel pipets suitable for 96- or 384-well format  
Barrier sterile low retention filter tips, 200 µl (GeneMate, cat. no. P-1237-200 for 96-well or CAPP, cat. no. 5030006C for 384-well format)  
Deep well 96-well plates (Eppendorf, cat. no. 951031909) or 384-well plates (Axygen, cat. no. P-384-240SQ-C-S)  
ELMI DTS-4 shakers (or equivalent shaking incubators capable of holding 96- or 384-well plates, shaking at 200 to 350 rpm, 37°C)  
Plate reader with optical density capabilities (600 nm) compatible with 96- or 384-well plate formats (e.g., Tecan Infinite M1000 Pro or BioTek Epoch 2)  
Vacuum aspirator setup (to which 200 µl pipet tips can be attached)

1. Following cell density determination for overnight cultures, add cells to wells:
  - a. For 96-well plate assays: At a final  $OD_{600} = 0.5$  (or equivalent to  $\sim 1 \times 10^7$  cells/ml) in 200 µl RPMI-1640 medium.
  - b. For 384-well plate assays: Add 1 µl overnight culture in 90 µl RPMI-1640 medium ( $\sim 2 \times 10^6$  cells per 90 µl).

*If starting a large number of wells with the same strain and medium, consider diluting cells into medium at the desired starting density and pipetting this mixture into wells from a sterile reservoir.*

2. Seal plates with Breathe-Easy® sealing membranes.

*Sealing plates reduces evaporation and prevents cross-contamination between wells.*


3. Shake plate at 37°C, 90 min at 250 rpm (96-well) or 350 rpm (384-well) in an ELMI incubator.
4. Remove membrane and aspirate medium.

*Change pipet tips on aspirator between wells of different strains and/or conditions. Take care not to scrape the bottom of the well with the tip when aspirating.*

5. Wash wells with 200 µl (96-well) or 50 µl (384-well) PBS. Aspirate PBS.
6. Add 200 µl (96-well) or 90 µl (384-well) fresh medium to each well.
7. Reseal plate with a new sealing membrane and shake at 37°C, 24 hr at 250 rpm (96-well) or 37°C, 24 hr at 350 rpm (384-well).
8. Remove membrane and aspirate medium from wells (Figure 2).

*Use clean pipet tips on the aspirator between wells of different strains and/or conditions. Take care not to scrape the bottom of the well with the tip when aspirating. Note any wells where biofilms have detached from the surface during aspiration so that they can be excluded from data analysis.*

9. Measure  $OD_{600}$  on a 96- or 384-well compatible plate reader.



**Figure 2** Typical *C. albicans* wild-type (WT) biofilm and negative control for the standard optical density assay in a 6-well plate format (modification of Basic Protocol 2). Shown is a 24 hr biofilm grown in Spider medium; the medium was aspirated immediately prior to visualization. Typical WT (left) and blank negative control (right) wells are shown. Placing the plate on a textured black background is recommended for the visualization of biofilms.

*The number of reads per well may vary based on the instrument. We recommend obtaining the average density of reads at five independent locations in each well in a 96-well plate or one read from the center of a 384-well plate. Be sure that the wells are still moist for the reading. Dry wells give inaccurate readings.*

10. Normalize data by subtracting OD<sub>600</sub> reading of an average blank well (containing medium alone) from each experimental and control well.

*The blank-subtracted OD<sub>600</sub> value of each experimental well (normally eight per condition in the 384-well format and six per condition in the 96-well format) is then normalized to the mean blank-subtracted OD<sub>600</sub> for the relevant control wells.*

11. Calculate mean, standard deviation, and statistical analyses (normally Student's unpaired two-tailed *t*-test assuming unequal variance) for each normalized data set.

### **2,3-BIS-(2-METHOXY-4-NITRO-5-SULFOPHENYL)-2H-TETRAZOLIUM-5-CARBOXANILIDE (XTT) REDUCTION ASSAY**

### **ALTERNATE PROTOCOL 1**

The 2,3-bis-(2-methoxy-4-nitro-5-sulfophenyl)-2H-tetrazolium-5-carboxanilide (XTT) assay is a colorimetric assay that detects metabolic activity by measuring the reduction of the tetrazolium salt reagent XTT; alternative versions of this assay use the tetrazolium salt reagent methylthiazolyldiphenyl-tetrazolium bromide (MTT; Krom et al., 2009; Krom, Cohen, McElhaney Feser, & Cihlar, 2007; Nett, Cain, Crawford, & Andes, 2011; Ramage, Vande Walle, Wickes, & Lopez-Ribot, 2001). This assay can be performed at the end of the standard optical density assay (Basic Protocol 2) and be combined with the inhibition and disruption optical density assays (Alternate Protocols 3 and 4). One or more variables may need to be optimized depending on the strains or species being used (e.g., exposure times, XTT concentrations). Note that the signal may not scale linearly between different strains and/or species, and much care should be taken when making these comparisons (Kuhn, Balkis, Chandra, Mukherjee, & Ghannoum, 2003). A preliminary time course assay (e.g., reading optical densities every 5 min for 1 hr after step 5 of this procedure) may be helpful to determine the exposure time best suited to a given strain.

#### **Additional Materials (also see Basic Protocol 2)**

- 0.5 mg/ml 2,3-bis-(2-methoxy-4-nitro-5-sulfophenyl)-2H-tetrazolium-5-carboxanilide (XTT; sodium salt, Sigma-Aldrich, cat. no. X4251) in PBS (make fresh prior to use and protect from light)



0.32 mg/ml phenazine methosulfate (PMS; Sigma-Aldrich, cat. no. P9625) in water (make fresh prior to use and protect from light)

Plate reader with optical density capabilities (492 nm) compatible with 96- or 384-well plate formats (e.g., Tecan Infinite M1000 Pro or BioTek Epoch 2)

1. Perform steps 1 to 7, Basic Protocol 2.
2. Approximately 15 to 30 min before the end of the 24-hr growth step, prepare fresh XTT and PMS solutions. Centrifuge before use to remove any insoluble materials and transfer carefully to a new tube. Protect from light.
3. Mix XTT and PMS at a 9:1 XTT/PMS ratio. Protect solution from light.
4. Remove membrane from plate, aspirate medium from wells.

*Consider measuring the optical density of the wells, as per step 9, Basic Protocol 2, as a second form of measurement before proceeding to step 5 of this protocol.*

5. Quickly add 100  $\mu$ l (96-well format) or 50  $\mu$ l (384-well format) of XTT/PMS mixture to wells while taking care to avoid disrupting the biofilm.

*Some protocols include one or more PBS washes prior to addition of the XTT/PMS solution. Based on our experience, omitting this wash step does not significantly affect results and reduces the likelihood of physical disruption of the biofilms. Should wash steps be performed, carefully track wells that have been physically disrupted so that they can be omitted from further analyses. This is a colorimetric assay based on the reduction of XTT; as such it is important to minimize the time between adding the XTT/PMS solution to the first and last wells on a given plate as large differences in reaction times between wells may skew results.*

6. Incubate plate for 30 min at 37°C in the dark.

*The incubation time may need to be increased or decreased depending on the strain(s) and conditions.*

7. Measure OD<sub>492</sub> on a 96- or 384-well compatible plate reader.

*Some protocols suggest transferring the liquid to a new plate prior to optical density measurement. Based on our experiences, this transfer is not necessary.*

## **ALTERNATE PROTOCOL 2**

### **CRYSTAL VIOLET ASSAY**

The crystal violet assay is another colorimetric assay that measures the uptake of a dye (crystal violet) by the biofilm. The amount of dye that ends up in the destaining solution serves as a proxy for the number of viable cells. This assay is less commonly used than the XTT reduction assay and requires more time and effort than either the XTT reduction assay or the standard optical density assay. The crystal violet assay can be performed at the end of the standard optical density assay and can be combined with the inhibition and disruption optical density assays (Jin, Yip, Samaranayake, Yau, & Samaranayake, 2003).

#### *Additional Materials (also see Basic Protocol 2)*

0.4% aqueous crystal violet solution (Sigma-Aldrich, cat. no. HT90132-1; store at room temperature for 4 to 6 months and protect from light)

95% ethanol

Plate reader with optical density capabilities (595 nm) compatible with 96- or 384-well plate formats (e.g., Tecan Infinite M1000 Pro or BioTek Epoch 2)

1. Perform steps 1 to 7, Basic Protocol 2.

2. After 24 hr growth step, remove membrane and aspirate medium from wells.
3. Wash wells twice with 200  $\mu$ l (96-well format) or 90  $\mu$ l (384-well format) PBS.

*Change pipet tips on aspirator between wells of different strains and/or conditions. Take care not to scrape the bottom of the well with the tip when aspirating or to physically disrupt the biofilm when adding liquid to the wells.*
4. Aspirate PBS and allow wells to air dry 45 min.
5. Stain with 110  $\mu$ l (96-well format) or 50  $\mu$ l (384-well format) 0.4% aqueous crystal violet solution 45 min.
6. Wash wells four times with 200  $\mu$ l (96-well format) or 90  $\mu$ l (384-well format) water.

*Change pipet tips on aspirator between wells of different strains and/or conditions. Take care not to scrape the bottom of the well with the tip when aspirating or to physically disrupt the biofilm when adding liquid to the wells.*
7. Destain wells 45 min in 200  $\mu$ l (96-well format) or 90  $\mu$ l (384-well format) 95% ethanol.
8. Transfer 100  $\mu$ l (96-well format) or 45  $\mu$ l (384-well format) 95% ethanol destaining solution from each well to a new plate.
9. Measure OD<sub>595</sub> on a 96- or 384-well compatible plate reader.

## INHIBITION OPTICAL DENSITY ASSAYS

The inhibition optical density assays modify the standard optical density assay to assess the ability of compounds to inhibit biofilm formation during the adherence and/or growth steps of the biofilm assay (Lohse et al., 2017). Solutions containing antifungal agents are included during the 90-min adherence step (adherence inhibition and sustained inhibition assays) and/or the 24-hr growth step (developmental inhibition and sustained inhibition assays). If it is feasible to perform only one version of the inhibition assays, we recommend that the sustained inhibition assay be chosen as it is the most likely to detect an effect. However, different compounds may have different effects in the three inhibition assays and as such we recommend that all three versions of the assay be performed when practical. Antifungal compounds are typically solubilized in DMSO. When preparing the antifungal solution in medium for the experiment, do not allow total DMSO concentration to exceed 2%. Prepare wells with DMSO at the same concentration as used in the assay, but without the compound of interest, to control for the effect of DMSO on the biofilm (this would serve as the control for normalization).

### *Additional Materials (also see Basic Protocol 2)*

Compound(s) to be tested  
Solvent(s), e.g., DMSO (Sigma-Aldrich, cat. no. D2650)

### *Adherence inhibition optical density assay*

- 1a. Include compound(s) being tested in medium when adding cells to medium during step 1, Basic Protocol 2.

*Compound is present for the 90-min adherence step, normal medium is present for the 24-hr growth step.*

- 2a. Perform steps 2 to 11, Basic Protocol 2.

### *Developmental inhibition optical density assay*

- 1b. Perform steps 1 to 5, Basic Protocol 2.

**ALTERNATE  
PROTOCOL 3**

Gulati et al.

**9 of 45**

- 2b. Include compound(s) being tested in medium when adding to wells during step 6, Basic Protocol 2.

*Normal medium is present for the 90-min adherence step, the compound is present for the 24-hr growth step.*

- 3b. Perform steps 7 to 11, Basic Protocol 2.

#### ***Sustained inhibition optical density assay***

- 1c. Include compound(s) being tested in medium when adding cells to medium during step 1, Basic Protocol 2.

*Compound is present for the 90-min adherence step and the 24-hr growth step.*

- 2c. Perform steps 2 to 5, Basic Protocol 2.

- 3c. Include compound(s) being tested in medium when adding to wells during step 6, Basic Protocol 2.

- 4c. Perform steps 7 to 11, Basic Protocol 2.

#### **ALTERNATE PROTOCOL 4**

#### **DISRUPTION OPTICAL DENSITY ASSAY**

The disruption optical density assay modifies the standard optical density assay to assess the ability of compounds of interest to disrupt an established, mature biofilm (Lohse et al., 2017). Although this protocol involves disrupting a 24-hr-old biofilm by exposing it to a compound for 24 hr, the age of the biofilm and the length of exposure to compounds can be modified as needed. This assay can be performed in 96- and 384-well formats, and we recommend the 384-well format whenever possible due to the reduced instances of physical (as opposed to chemical) disruption of biofilms during the disruption step. This protocol involves the removal of the existing medium from the biofilms and the addition of fresh medium containing the compound(s) being tested. When screening large compound libraries, it may be worth considering adding the compound(s) directly to the existing medium rather than removing the existing medium and adding fresh medium with the compound(s). Antifungal compounds are typically solubilized in DMSO. When preparing the antifungal solution in medium for the experiment, do not allow the final DMSO concentration to exceed 2%. Prepare wells with DMSO, at the same concentration as used in the assay but without the compound of interest to control for the effect of DMSO on the biofilm (this would serve as the control for normalization).

#### ***Additional Materials (also see Basic Protocol 2)***

Compound(s) to be tested


Solvent(s), e.g., DMSO (Sigma-Aldrich, cat. no. D2650)

1. Perform steps 1 through 7, Basic Protocol 2.

*We normally use 8 wells per compound for this assay. Since the same strain is normally used for all compounds being tested, the density of the overnight culture is not determined and instead 1  $\mu$ l overnight culture (equivalent to  $\sim 2 \times 10^6$  cells/ml) is added to the medium in each well in the 384-well version of this assay.*

2. After the 24-hr growth step, remove membrane and carefully aspirate medium in groups of 6 to 12 wells.

*Change pipet tips on aspirator between wells of different strains and/or conditions. Aspirate by slowly lowering tip into well until all liquid has been removed; take care not to scrape the bottom of the well with tip when aspirating.*



**Figure 3** Standard aspiration techniques and locations in microtiter plates. **(A)** When aspirating, place the tip against the wall of the well above the surface of the liquid (far left). Slowly lower the tip, taking care to keep it against the well wall (middle left). Stop aspirating when all liquid is gone, which should occur before the tip reaches the bottom of the well and is just above the biofilm (middle right). Do not lower the tip all the way to the bottom of the well as this will physically disrupt any biofilm present (far right). **(B-C)** Locations to use when aspirating and adding liquid to a well from a 96-well (B) and 384-well (C) microtiter plate. When aspirating at position “1”, add liquid back at position “2”. We normally aspirate from one of the corners of a square well of a 384-well microtiter plate.

3. Add 200  $\mu\text{l}$  (96-well) or 90  $\mu\text{l}$  (384-well) medium containing the compound of interest to the wells. Add medium slowly to the side of the well opposite the side from which medium was aspirated (Figure 3).

*Add medium slowly to reduce physical disruption of the biofilm, especially in 96-well format. Aspirate and add back medium to one group of 6 to 12 wells at a time in order to avoid exposing biofilms to air and/or desiccation for extended periods of time. Note any wells where the biofilm was disrupted during the aspiration/addition process.*

4. Reseal plate and shake at 37°C an additional 24 hr at 250 rpm (96-well) or 350 rpm (384-well).
5. After 24 hr, read plate and analyze as per steps 8 to 11, Basic Protocol 2.

## DISPERSAL OPTICAL DENSITY ASSAYS

The dispersal optical density assays modify the standard optical density assay to measure the cells that are dispersed during biofilm growth (Lohse et al., 2017; Nobile et al., 2014). These assays utilize optical density to measure the dispersed cells at 24, 48, and 60 hr time points. There are two assays that can be used to assess cell dispersal during biofilm formation, the sustained dispersal assay and the standard dispersal assay. The sustained dispersal assay measurements are obtained while the biofilm grows in the original medium with no additional medium added, requiring a separate set of samples for each time point. In the standard dispersal assay, in contrast, fresh medium is added back to the biofilm after each measurement. These dispersal assays can be customized to study biofilm dispersal differences between different *C. albicans* strains or to measure the effect of a compound on dispersal during biofilm formation. The assays can be performed in 96- and 384-well plate formats. We recommend the 96-well format whenever possible due to the ease of visually observing dispersed cells floating above the biofilm in the wells. Although slightly more involved than the standard optical density assay, from which it is derived, it is still possible to screen large numbers of mutant strains or compound libraries for defects in biofilm dispersal using these assays. We recommend six wells per strain or condition in 96-well plate format; this allows for evaluation of sixteen conditions per 96-well plate (including blanks and controls).

## Materials

See Basic Protocol 2

### ***Sustained dispersal optical density assay***

- 1a. Perform steps 1 through 7, Basic Protocol 2.

*A separate set of replicates is required for each time point, as once the dispersed cells are measured, medium cannot be added back to the well containing the biofilm.*

- 2a. After the 24-hr growth step, remove membrane and carefully remove medium with the dispersed cells from the first set of replicates and add this medium to a clean 96- or 384-well plate.

*Take care not to scrape the bottom of the well or to disturb the biofilm on the walls of the well with the tip during this process.*

- 3a. Homogenize the removed medium using a multichannel pipet.
- 4a. Measure OD<sub>600</sub> of removed medium (containing dispersed cells) and measure OD<sub>600</sub> of the biofilm, from which medium was removed, on a 96- or 384-well compatible plate reader.
- 5a. Repeat steps 2a to 4a for the second set of replicates at the 48-hr time point, followed by the third set of replicates at the 60-hr time point.
- 6a. Normalize data at each time point by subtracting the OD<sub>600</sub> reading of an average blank well (containing medium alone) from each experimental and control well from that time point.
- 7a. Divide blank-subtracted OD<sub>600</sub> value for medium from each experimental well by the blank-subtracted OD<sub>600</sub> value of biofilm of the correlated well to obtain normalized OD<sub>600</sub> values.

*This provides a measurement of dispersed cells per biofilm and takes into account any variations in biofilms between wells.*

- 8a. Divide normalized OD<sub>600</sub> value for each of the experimental conditions by the normalized OD<sub>600</sub> value for control wells. Perform step 13 below.

### ***Standard dispersal optical density assay***

- 1b. Perform steps 1 through 7, Basic Protocol 2.

- 2b. After the 24-hr growth step, remove membrane and carefully remove medium with dispersed cells from wells with biofilms and add this medium to a clean 96- or 384-well plate.

*Take care not to scrape the bottom of the well or disturb biofilm on the walls of the well with tip during this process.*

- 3b. Add 200  $\mu$ l (96-well) or 90  $\mu$ l (384-well) fresh medium to the original wells. Add medium slowly to the side of the well opposite the side from which medium was removed (Figure 3).

*Add medium slowly to reduce physical disruption of the biofilm, especially in 96-well format.*

- 4b. Reseal plate and shake at 37°C, 250 rpm (96-well) or 350 rpm (384-well).
- 5b. Homogenize removed medium using a multichannel pipet.
- 6b. Measure OD<sub>600</sub> of removed medium (containing dispersed cells) on a 96- or 384-well compatible plate reader.

- 7b. After an additional 24 hr, repeat steps 2b to 6b.
- 8b. After an additional 12 hr, remove membrane and carefully remove medium with dispersed cells from wells with biofilms and add medium to a clean 96- or 384-well plate.
- 9b. Homogenize removed medium using a multichannel pipet.
- 10b. Measure OD<sub>600</sub> of removed medium (containing dispersed cells) and measure OD<sub>600</sub> of biofilm well from which medium was removed using a 96- or 384-well compatible plate reader.
- 11b. Normalize data at each time point by subtracting the OD<sub>600</sub> reading of an average blank well (containing medium alone) from each experimental and control well from that time point.
- 12b. Divide blank-subtracted OD<sub>600</sub> value for medium from each experimental well by the blank-subtracted OD<sub>600</sub> value of biofilm of the correlated well to obtain normalized OD<sub>600</sub> values.

*This provides a measurement of dispersed cells per biofilm at a given time point and takes into account any variations in biofilms between wells.*
13. Calculate mean, standard deviation, and statistical analyses (typically use a Student's unpaired two-tailed *t*-test assuming unequal variance) for each normalized data set.

#### **CO-CULTURING AND ANALYZING *CANDIDA* MIXED-SPECIES BIOFILMS USING THE STANDARD OPTICAL DENSITY ASSAY**

#### **BASIC PROTOCOL 4**

This protocol is optimized for culturing *C. albicans* mixed-species biofilms with *Escherichia coli*, however, the culture conditions indicated can be applied or adapted for growing other microbes together with *C. albicans*. This protocol is a refined version of the protocol reported by Fox et al. (2104), for co-culturing *C. albicans* with different anaerobic and aerobic bacteria. In the following protocol, *E. coli* was chosen for co-culture with *C. albicans* as both species are already known to interact in the context of a biofilm (Fox et al., 2014), are prevalent in the human gut as well as other parts of the body, and can be easily cultured in the laboratory. The 6- and 12-well polystyrene plate biofilm assay described in this protocol can be modified to use 96- and 384-well plate formats for higher throughput screening. Likewise, the protocol can be modified to accommodate the requirements of many of the other protocols derived from the standard optical density assay. We note that this protocol does not assess the relative abundance of each of the two species (in terms of cell numbers) within the mature biofilm or the viability of either species within that biofilm. If these determinations are of interest, we recommend combining this assay with one or both versions of the cell enumeration assay.

##### ***Additional Materials* (also see *Basic Protocol 2*)**

- E. coli* or other microbial strain(s) of interest
- Luria-Bertani (LB) medium for *E. coli* growth (or appropriate medium for other microbes under study)
- Bacto heart infusion (BHI) medium (BD, manufacturer no. B237500), supplemented with 5% FBS (BHI-FBS)
- Plates for *E. coli* growth
- 6-well or 12-well microtiter plates, non-tissue culture-treated (Falcon brand, manufacturer no. 351146 and 351143)
- 5-ml serological pipets (VWR, cat. no. 89130-896)

Gulati et al.

---

**13 of 45**

1. Streak *E. coli* strains onto LB plates from glycerol stocks 2 to 5 days in advance of assay. Incubate plates at 37°C, 24 hr.

*Plates containing E. coli can be stored at 4°C for up to 1 week.*

**For small numbers of strains/large volumes**

- 2a. Inoculate a single *E. coli* colony in 25 ml LB medium in a 250-ml flask from a 2 to 5 day old plate. Grow at 37°C with shaking overnight.

**For large numbers of strains/small volumes**

- 2b. Inoculate a single *E. coli* colony in 4 ml LB medium in a 20-ml test tube from a 2 to 5 day old plate. Grow at 37°C with shaking or using a roller drum overnight.
3. After 12 to 16 hr, dilute overnight culture 1:20 or 1:40 and measure OD<sub>600</sub> using cuvettes and a spectrophotometer or 96-well plates and a plate reader.

*If assays will be set up more than 16 hr after overnight cultures were started, remove overnight cultures from the shaker or roller drum in the morning and allow to sit at room temperature until time of use (cultures should not sit for more than 2 hr at room temperature). Vortex cultures immediately before use. However, if possible, it is better to time the experiment such that the strains will be used prior to 16 hr of overnight growth.*

4. For culturing *C. albicans* strains, follow Basic Protocol 1.
5. Following cell density determination for overnight cultures, add cells to wells at a final OD<sub>600</sub> = 0.5 (*C. albicans*) or OD<sub>600</sub> = 0.0125 (*E. coli*); equivalent to ~1 × 10<sup>7</sup> cells/ml for each species in 4 ml (for 6-well plate assays) or 2 ml (for 12-well plate assays) BHI-FBS medium.

*This protocol seeds the same number of bacterial and fungal cells together at the start of the experiment in a 1:1 ratio, however this ratio can be altered as needed. When working with a species for the first time, we recommend the use of a serial dilution and colony forming unit (CFU) count to correlate OD<sub>600</sub> values with an approximate number of cells/ml so that the 1:1 ratio can be achieved (see Alternate Protocol 5). It is also recommended that each species be grown as biofilms individually in BHI-FBS as a single-species biofilm control to compare to the mixed-species biofilms.*

6. Shake plate at 37°C, 90 min at 200 rpm (6-well plate) or 250 rpm (12-well plate) in an ELMI incubator.
7. Aspirate medium and wash wells with 2 ml (12-well) or 4 ml (6-well) PBS.
8. Aspirate PBS and add 2 ml (12-well) or 4 ml (6-well) of fresh medium to each well.
9. Shake plate at 37°C, 24 hr at 200 rpm (6-well) or 250 rpm (12-well).
10. Aspirate medium from wells.

*Change pipet tips on aspirator between wells of different strains and/or conditions. Take care not to scrape the bottom of the well with the tip when aspirating. Note any wells where biofilms detached from the surface during aspiration so that they can be excluded from data analysis.*

11. Measure OD<sub>600</sub> on a 6- or 12-well compatible plate reader.

*The number of reads per well may vary based on the instrument. We recommend taking the average density of five reads in each section of a five-by-five grid (i.e., twenty-five independent locations) in each well in a 6-well plate or 12-well plate. Be sure that the wells are still moist for the reading as dry wells give inaccurate readings.*

12. Normalize data by subtracting OD<sub>600</sub> reading of an average blank well (containing medium alone) from each experimental and control well. The blank-subtracted OD<sub>600</sub> value of each experimental well is then normalized to the mean blank-subtracted OD<sub>600</sub> for relevant control wells.
13. Calculate standard deviation and statistical analyses (typically use a Student's unpaired two-tailed *t*-test assuming unequal variance) for each normalized data set.

### COLONY-FORMING UNITS ASSAY

The starting cell concentration for a biofilm assay as well as the ratio of fungal or bacterial species used to seed mixed-species biofilms significantly affects the nature of the resulting biofilms. As such, accurate quantitation of microbial samples is an essential precursor to growing different species in a biofilm. The CFU assay presented here is a basic method to correlate an unknown OD<sub>600</sub> concentration with a specific number of viable cells.

#### *Materials*

*C. albicans*, *E. coli*, or other microbial strains of interest

Luria-Bertani (LB) medium (or appropriate medium for other microbes under study)

Dulbecco's phosphate-buffered saline (PBS; calcium and magnesium salt free), sterile filtered

YEPD agar plates (2% Bacto™ peptone, 2% dextrose, 1% yeast extract, 2% agar; store at room temperature)

Plates for *E. coli* growth

1.5-ml microcentrifuge tubes

1. Culture strains following steps 1 through 3, Basic Protocol 1 (for fungi) or Basic Protocol 4 (for bacteria) with modifications as appropriate for the species in question.
2. Following cell density determination for overnight cultures, dilute cells to a final OD<sub>600</sub> = 0.5 (for fungi) or OD<sub>600</sub> = 0.0125 (for *E. coli*); equivalent to ~1 × 10<sup>7</sup> cells/ml for each species in 1 ml PBS.

*These OD<sub>600</sub> values are an approximate estimate for a cell concentration of 1 × 10<sup>7</sup> cells/ml to allow for a serial dilution series that will more accurately correlate OD<sub>600</sub> with viable cells per ml.*

3. Perform a series of six sequential ten times dilutions by adding 100 μl cell-containing solution to 900 μl PBS. Vortex thoroughly and change pipet tips after each dilution.

*The resulting dilution series should range from a high of 1 × 10<sup>7</sup> to a low of 10 cells/ml.*

4. Vortex tubes thoroughly and plate 100 μl from each of tubes 3 through 7 on YEPD plates.

*We recommend performing at least two separate dilution series, each based on separate initial OD<sub>600</sub> measurements, for each strain.*

5. Incubate plates at 30°C, 24 to 72 hr or other conditions as appropriate for the strain(s).
6. Count colonies on the plate and use the serial dilution factor to estimate number of cells in the original culture for each replicate in a set.



*When counting colonies, plates should be examined within 24 hr to ensure that colonies can be distinguished before they overgrow if too many colonies were plated. Plates should also be examined after 48 to 72 hr to allow for the scoring of slower growing strains or species. This set of time points allows for adequate growth while ensuring that colonies are easily distinguished from one another.*

**BASIC  
PROTOCOL 5**

**SILICONE SQUARE ASSAY**

For some protocols, it may be of interest to form *C. albicans* biofilms on silicone squares rather than (or in addition to) polystyrene plates to test the effects of an alternate surface on biofilm formation (Nobile, Andes, et al., 2006; Nobile et al., 2012; Nobile & Mitchell, 2005; Nobile, Nett, et al., 2006). Silicone is a common material used for implanted medical devices, such as catheters, heart pumps, and surgical reconstructive components. Biofilms formed using the silicone square assay can be imaged using confocal microscopy or weighed to determine the dry weight (biomass) of the biofilm formed. For the latter, it is essential to pre-weigh the silicone squares prior to use.

**Materials**

RPMI-1640 medium with L-glutamine and 3-(N-morpholino)propanesulfonic acid (MOPS), without sodium bicarbonate, pH 7.0 (see recipe) or other medium of interest

Dulbecco's phosphate-buffered saline (PBS; calcium and magnesium salt free), sterile filtered

200 proof ethanol (CAS 64-17-5; e.g., Thermo Fisher Scientific, cat. no. BP2818100)

Autoclaved silicone squares (1.5 × 1.5 mm) cut from medical-grade silicone sheets (Bentec Medical, cat. no. PR72034-06N)

12-well microtiter plates, non-tissue culture-treated (Falcon brand, manufacturer no. 351143)

Dissecting forceps with a fine tip (VWR, cat. no. 82027-404)

ELMI DTS-4 shakers (or equivalent shaking incubators capable of holding 12-well plates, shaking at 200 rpm, 37°C)

1. Place one autoclaved silicone square per well into a 12-well plate using sterile forceps.


*Weigh each silicone square before placing them in the wells if you plan to use them for dry weight assays. Track which square went into each well. Include positive (e.g., wild-type strain) and negative (e.g., known biofilm-defective mutant strain) control wells in addition to a blank (uninoculated) control well.*

2. Following cell density determination for overnight cultures, add cells to silicone squares at a final OD<sub>600</sub> = 0.5 (equivalent to ~1 × 10<sup>7</sup> cells/ml) in 2 ml RPMI-1640 medium.

*Care must be taken to note the direction (facing up) of the silicone square in the well, as this is the surface where cells will start to adhere. In all future steps, this side should always face upwards.*

3. Shake plate at 37°C, 90 min at 200 rpm in an ELMI incubator.
4. Add 2 ml PBS to wells in a new 12-well plate. Gently lift silicone square from the edges using forceps and place in well containing PBS.

*Carefully lift silicone squares by holding just the edges, so as to not disturb the adhered cells.*



**Figure 4** Typical *C. albicans* wild-type (WT) biofilm and negative control for the silicone square assay (Basic Protocol 5). The biofilm was allowed to develop on a silicone square for 24 hr in Spider medium in a 6-well plate. The biofilm is clearly visible on the surface of the silicone square with WT (left). A blank silicone square is shown as a negative control (right).

5. Wash silicone square in PBS by holding the edges using forceps and by gently agitating the square by moving it up and down in the PBS solution two to three times.
6. Place silicone square into the well of a new 12-well plate containing 2 ml fresh RPMI-1640 medium.
7. Shake plate at 37°C, 24 hr at 200 rpm in an ELMI incubator (Figure 4).

*The resulting biofilm (Figure 4) can be visualized using confocal microscopy or used to measure dry weight (biomass) of the biofilm. To measure biofilm dry weights after growth on silicone squares, simply weigh the silicone square after allowing the squares containing the biofilm to dry in an ELMI at 37°C, 4 hr; no shaking is necessary for this step. Subtract the weight of the silicone square without the biofilm from the weight of the silicone square containing the biofilm to obtain the biofilm dry weight.*

## CELL ADHESION ASSAY

Initial cell adhesion to a surface is a crucial step in normal biofilm development and is the first step of biofilm formation. The cell adhesion assay measures the number of cells that adhere to an abiotic surface, such as a microtiter plate, to seed a biofilm. This assay mimics the adherence step of the standard optical density assay and can be used to measure the effect of a compound(s) on cell adhesion or compare the ability of strains to adhere to a surface (Lohse et al., 2017; Winter et al., 2016). Note that this assay has reduced throughput compared to the standard optical density assay as it involves serial dilutions and plating to count CFUs to determine the number of adhered cells. We recommend that this assay be performed in a 96-well plate format, as this format provides an ideal surface area for cell adherence. We recommend the use of four wells per strain or condition in the 96-well plate format, allowing for the evaluation of 24 conditions per 96-well plate (including blanks and controls). The assay can be customized for 6-well and 12-well plate formats, however the serial dilutions should be optimized to correlate with the well surface area of the 6- and 12-well plates.

### *Additional Materials* (also see Basic Protocol 2)

Compound(s) to be tested

YEPD agar plates (2% Bacto™ peptone, 2% dextrose, 1% yeast extract, 2% agar; store at room temperature)

1.5-ml microcentrifuge tubes

1. Perform steps 1 through 4, Basic Protocol 2 for the 96-well plate.

2. Wash wells twice with 200  $\mu$ l PBS. Aspirate PBS.  
*This step removes non- and weakly adhered cells.*
3. Add 200  $\mu$ l PBS and vigorously resuspend adhered cells.  
*Scrape the bottoms and edges of each well with a pipet tip to dislodge cells that remain adhered to the surface. Use the pipet to homogenize the suspension.*
4. Perform serial dilutions in PBS. Make first dilution by adding 200  $\mu$ l PBS with resuspended cells to 1800  $\mu$ l PBS to obtain a ten times dilution. Vortex solution thoroughly.
5. Perform a second dilution by adding 100  $\mu$ l of the solution from step 4 to 900  $\mu$ l PBS. Vortex solution thoroughly and plate 100  $\mu$ l solution onto a YEPD plate.  
*A 1000-fold serial dilution normally provides an optimal number of cells for counting for wild-type conditions. The serial dilution may need to be optimized for different compound(s) and/or strains, depending on their effect on cell adherence. If cell adherence is too high, such that single colonies cannot be effectively counted, perform another 10 $\times$  dilution. If cell adherence is too low, such that not enough colonies are obtained, plate 100  $\mu$ l from the first dilution tube.*
6. Incubate YEPD plates at 30°C, 24 to 72 hr.
7. Count colonies on each plate and use the serial dilution factor to determine the number of adherent cells for each replicate in a set.  
*When counting colonies, plates should be examined within 24 hr to ensure that colonies can be distinguished before they overgrow if too many cells were plated. Plates should also be examined after 48 to 72 hr to allow for the scoring of slow growing mutant strains. This set of time points allows colonies adequate time to grow while ensuring that they are easily distinguished from one another.*
8. Calculate mean, standard deviation, and statistical analyses (typically use a Student's unpaired two-tailed *t*-test assuming unequal variance) for each data set.

**BASIC  
PROTOCOL 7**

**DRY WEIGHT ASSAY**

This assay measures the biomass or dry weight of *C. albicans* biofilms grown on the bottom of polystyrene plates in 6-well and 12-well formats (Hawser & Douglas, 1994; Nobile, Andes, et al., 2006; Nobile et al., 2012; Nobile, Nett, et al., 2006). We typically use three replicates per strain or condition. We do not recommend the 96-well or 384-well formats for this assay, as the total mass of biofilm harvested from wells from these plates is much less than for 6- or 12-well plates, and the reduced starting biomass makes it difficult to assess the effects of treatment with compound(s) or to assess differences between distinct strains (e.g., mutant strains of interest). Although this assay is effective in detecting severe biofilm defects, the dry weight assay is less sensitive relative to other more recently developed biofilm assays, such as the standard optical density assay, and thus is not ideal for detecting minor differences in biofilm formation. The dry weight assay is also notoriously ineffective at detecting enhanced biofilm formation and should not be used to assess such strains. This assay also has lower throughput than many of the other biofilm assays presented. For these reasons, we recommend the dry weight assay as a secondary assay to further assess candidates identified in other *in vitro* biofilm assays.

**Materials**

RPMI-1640 medium with L-glutamine and 3-(*N*-morpholino)propanesulfonic acid (MOPS), without sodium bicarbonate, pH 7.0 (see recipe)  
Dulbecco's phosphate-buffered saline (PBS; calcium and magnesium salt free), sterile filtered

70% ethanol

Mixed cellulose esters membrane (Millipore, cat. no. AAWG02500)

Millipore filter device for 25-mm discs (e.g., Millipore, cat. no. XX10025400 or equivalent)

6-well or 12-well microtiter plates, non-tissue culture-treated (Falcon brand, manufacturer no. 351146 and 351143)

Vacuum aspirator setup (to which 1000  $\mu$ l tips can be attached)

Dissecting forceps with a fine tip (VWR, cat. no. 82027-404)

Analytical scale

ELMI DTS-4 shakers (or equivalent shaking incubators capable of holding 6- and 12-well plates, shaking at 200 to 250 rpm, 37°C)

1. Following cell density determination for overnight cultures, add cells to wells at a final  $OD_{600} = 0.5$  (or equivalent to  $\sim 1 \times 10^7$  cells/ml) in 2 ml (for 12-well plate assays) or 4 ml (for 6-well plate assays) RPMI-1640 medium.

*Although the protocol listed here uses RPMI-1640 medium, other biofilm inducing media can alternatively be used. The authors recommend leaving one well blank (containing medium alone with no cells) as a contamination control and having at least three replicates per condition.*

2. Shake plate at 37°C, 90 min at 200 rpm (6-well plate) or 250 rpm (12-well plate) in an ELMI incubator.
3. Aspirate medium and wash wells with 2 ml (12-well) or 4 ml (6-well) PBS.


*Change pipet tips on aspirator between wells of different strains and/or conditions. Take care not to scrape the bottom of the well with tip when aspirating.*

4. Aspirate PBS and add 2 ml (12-well) or 4 ml (6-well) fresh medium to each well.
5. Shake plate at 37°C, 24 hr at 200 rpm (6-well) or 250 rpm (12-well).
6. Remove plate and aspirate medium from wells.
7. Vigorously disrupt biofilm and resuspend all cells in 2 ml (12-well) or 4 ml (6-well) PBS.

*Scrape the bottoms and edges of each well with a pipet tip to dislodge biofilm and cells that remain adhered to the surface. Use the pipet to homogenize the biofilm suspended in PBS.*

8. Clean Millipore filtration device (or equivalent) parts with 70% ethanol and rinse thoroughly with sterile water.
9. Assemble filtration device by first placing the plastic ring on the waste flask to maintain a tight seal. Then add, in this order: Wire filter sieve, mixed cellulose ester membrane, followed by the funnel. Lastly, attach clamp to keep all components together (Figure 5).
10. Attach vacuum aspirator setup to the filtration device.
11. Turn vacuum on and transfer disrupted biofilm of one well (from step 7) to the funnel to adhere to the membrane. Apply biofilm solution towards the center of the membrane.

*If additional biofilm is left behind in the well, an additional 1 ml PBS can be added to the well to suspend the remaining biofilm. Then transfer this additional disrupted biofilm to the membrane.*



**Figure 5** Assembly of the vacuum filtration apparatus used for the dry weight assay (Basic Protocol 7). The filtration device is assembled by placing the rubber cork on a flask (A), followed by the plastic O-ring (B), the wire filter sieve (C), and the mixed cellulose ester membrane (D; in this order from bottom to top). The funnel (E) is then added on top and the clamp (F) is attached to keep all components together. The clamp should cover the top of the glass piece (A) as well as the base of the funnel (E).

12. Turn off vacuum, disassemble filtration device, and carefully remove wet membrane using forceps. Only use the edges of the membrane to avoid disturbing the accumulated cells.

13. Place membrane in a clean well of a 6-well plate.

*Use a 6-well plate for this step as the membrane is too large to sit flat in a 12-well plate.*

14. Repeat steps 8 through 14 for all wells, including blank control wells.

*The filter paper from the blank wells will be used to determine the average weight of the mixed cellulose ester membranes containing any medium components. Note that the weights of the individual ester membranes are nearly identical and thus it is not necessary to pre-weigh each ester membrane prior to use.*

15. Place the 6-well plates with the cellulose ester membranes in an ELMI at 37°C, 4 hr; no shaking is necessary for this step.

*It is important to allow the membranes to dry thoroughly and uniformly as uneven wetness can lead to inaccurate weight measurements.*

16. Weigh cellulose ester membranes containing biofilm cells individually. Subtract the average of the blank wells from all replicates.

17. Calculate mean, standard deviation, and statistical analyses (typically use a Student's unpaired two-tailed *t*-test assuming unequal variance) for each data set.

## CELL ENUMERATION ASSAY

This assay determines the viability of cells within a biofilm after exposure to a compound of interest, such as an antifungal agent, to assess if the compound has microbicidal

properties against the cells within the biofilm. This protocol involves growing a wild-type biofilm and exposing the biofilm to the compound being tested for 24 hr, although the assay can be modified to use different exposure times based on the needs of the researcher. The addition of Lethen broth, which enhances biofilm homogenization, avoids the difficulties of obtaining inaccurate colony counts when biofilms are not effectively broken apart, which can occur when biofilms are grown in certain media (e.g., RPMI-1640). We recommend that this assay be performed in 96-well plate formats as this allows for a larger number of live cells present compared to biofilms formed in 384-well plate formats. We recommend four wells per strain or condition in 96-well plate format, allowing for evaluation of 24 conditions per 96-well plate (including blanks and controls). This assay can also be used to assess the relative abundance and viability of the individual species present in a mixed-species biofilm (Fox et al., 2014). When working with mixed-species biofilms, it is sometimes possible to get efficient homogenization without the use of Lethen broth depending on the medium and species used (this should be microscopically verified for each specific case). Should Lethen broth prove necessary, we note that it is effective at breaking apart both fungal and bacterial biofilms.

**Additional Materials** (also see *Basic Protocol 2* and *Alternate Protocol 4*)

Compound(s) to be tested

Solvent(s), e.g., DMSO (Sigma-Aldrich, cat. no. D2650)

Lethen broth (Difco brand, manufacturer no. 268110 or see recipe)

Sodium thiosulfate (Thermo Fisher Scientific, cat. no. S445–500; prepare fresh 10% solution and dilute as required)

1.5-ml microcentrifuge tubes

1. Perform steps 1 through 7, Basic Protocol 2 followed by steps 2 through 4, Alternate Protocol 4.
2. After 24 hr, carefully remove plate from the shaker and remove membrane.
3. Carefully aspirate medium in groups of 6 to 12 wells and gently add 200  $\mu$ l PBS to wash biofilms.

*At this point, the biofilms are very brittle. To wash the biofilms, add PBS to the wells in a drop-wise manner from one corner of the well to avoid mechanical disruption (Figure 3). Single-channel pipets are recommended over multi-channel pipets for better control at this step.*

4. Carefully aspirate PBS and add 200  $\mu$ l fresh PBS to the wells. Vigorously disrupt biofilm to resuspend cells.

*Scrape the bottoms and edges of each well with a pipet tip to dislodge the biofilm and cells that remain adhered to the surface. Use the pipet to homogenize the biofilm suspended in PBS.*

5. Dilute resuspended biofilm into 2 ml Lethen broth supplemented with 0.1% sodium thiosulfate and vigorously vortex to further homogenize the biofilm.
6. Perform serial dilutions by setting up a series of four tubes for ten times dilutions in PBS for each well. Add 100  $\mu$ l cells resuspended in Lethen broth and sodium thiosulfate to 900  $\mu$ l PBS to obtain a ten times dilution (tube 1). Vortex thoroughly.
7. Perform a second ten times dilution, by adding 100  $\mu$ l tube 1 to 900  $\mu$ l PBS (tube 2). Vortex thoroughly.
8. Repeat step 7 twice to obtain tubes 3 and 4 (diluting tube 2 into tube 3 and tube 3 into tube 4). Vortex tube 4 thoroughly and plate 100  $\mu$ l solution on a YEPD plate;

there should be one plate per replicate, thus a total of four plates for each strain or condition.

*These dilutions have been optimized for a wild-type *C. albicans* biofilm with no treatment (control). When testing a new compound or strain, it is recommended that 100  $\mu$ l from tube 3 should also be plated on a separate YEPD plate. This ensures that colonies will be obtained even if there is a higher instance of cell death relative to the wild-type strain.*

9. Incubate plates at 30°C, 24 to 72 hr.
10. Count colonies on the plate within 24 hr and between 48 to 72 hr and use the serial dilution factor to estimate the number of cells in the biofilm for each replicate in a set.

*When counting colonies, plates should be examined within 24 hr to ensure that colonies can be distinguished before they overgrow if too many colonies were plated. Plates should also be examined between 48 to 72 hr to allow for the scoring of slow growing strains. This set of time points allows for adequate growth while ensuring that colonies are easily distinguished from one another.*

11. Calculate mean, standard deviation, and statistical analyses (typically use a Student's unpaired two-tailed *t*-test assuming unequal variance) for each data set.

## **ALTERNATE PROTOCOL 6**

### **CELL ENUMERATION USING FLORESCENCE MICROSCOPY**

The survival of cells in *C. albicans* biofilms treated with antifungal agents can also be quantified using fluorescent live/dead stains rather than colony counts of serial dilutions. This protocol has been adapted from LaFleur, Kumamoto, and Lewis (2006) to work with the disruption variant of the standard optical density assay and requires less time and resources compared to the serial dilution method. Additionally, the visualization step better accounts for any cells that are still adhered to one another, which would be overlooked in the serial dilution version of this protocol. This protocol requires access to a microscope with fluorescent capabilities. We recommend that this assay be performed in 96-well plate format and that three wells per strain or condition be tested; this allows for evaluation of 32 conditions per 96-well plate (including blanks and controls).

***Additional Materials*** (also see *Basic Protocol 2, Basic Protocol 8, and Alternate Protocol 4*)

- 5 mg/ml fluorescein diacetate (Sigma-Aldrich, cat. no. F7378-5G) in DMSO (store at  $-20^{\circ}\text{C}$  per manufacturer's recommendations and protect from light)
- Microscope slides (e.g., Thermo Fisher Scientific, cat. no. 12-550-123 or equivalent)
- Coverslips (e.g., Thermo Fisher Scientific, cat. no. 12-542-B or equivalent)
- Fluorescent microscope with ability to detect green fluorescent protein (GFP)
- Microcentrifuge (Eppendorf miniSpin plus or equivalent)
- ImageJ software [National Institutes of Health (NIH)] or equivalent software

1. Perform steps 1 through 7, Basic Protocol 2.
2. Prepare RPMI-1640 with the antifungal agent of interest and 100  $\mu\text{g}/\text{ml}$  fluorescein diacetate (live/dead stain).
3. After the 24-hr growth step, remove membrane and carefully aspirate medium in groups of 6 to 12 wells.

*Change pipet tips on aspirator between wells of different strains and/or conditions. Aspirate by slowly lowering tip into well until all liquid has been removed. Take care not to scrape the bottom of the well with the tip when aspirating.*

4. Add 200  $\mu$ l medium containing the compound of interest and fluorescein diacetate to the wells. Add medium slowly to the side of the well opposite the side from which medium was aspirated (Figure 3).

*Add medium slowly to reduce physical disruption of the biofilm, especially in 96-well plate format. Aspirate and add back medium to one group of 6 to 12 wells at a time in order to avoid exposing biofilms to air and/or desiccation for extended periods of time. Note any wells where the biofilm was disrupted during the aspiration/addition process.*

5. Reseal plate and shake at 37°C, an additional 24 hr at 250 rpm.
6. After 24 hr, carefully remove plate from the shaker and remove membrane.
7. Carefully aspirate medium in groups of 6 to 12 wells and gently add 200  $\mu$ l PBS to wash biofilms (Figure 3).

*At this point, the biofilms are very brittle. To wash the biofilms, add PBS to these wells in a drop-wise manner from one corner of the well to avoid mechanical disruption. Single-channel pipets are recommended over multi-channel pipets for better control at this step.*

8. Carefully aspirate PBS and add 200  $\mu$ l fresh PBS to the wells. Vigorously disrupt biofilm to resuspend the cells.

*Scrape the bottoms and edges of each well with a pipet tip to dislodge the biofilm and cells that remain adhered to the surface. Use the pipet to homogenize the biofilm suspended in PBS.*

9. If Lethen broth is needed to break apart the biofilm, dilute resuspended biofilm into 2 ml Lethen broth supplemented with 0.1% sodium thiosulfate and vigorously vortex to further homogenize the biofilm. Otherwise, proceed to step 10.
10. Vortex vigorously and transfer disrupted biofilm from each well to a 1.5-ml microcentrifuge tube.
11. Centrifuge biofilm-containing tubes in a microcentrifuge 5 min at 9000  $\times$  g to separate cells from excess dye in the solution.
12. Discard supernatant and resuspend pellet in 200  $\mu$ l PBS.
13. Take 1  $\mu$ l resuspended pellet and place on a clean microscope slide. Place glass coverslip on top of the sample.

*Make one slide per replicate well, for a total of three slides for each strain or condition being tested.*


14. Visualize slide using a fluorescence microscope, utilizing the GFP settings.
15. Use 20 $\times$  objective lens to focus on the cells in the biofilm and take an image in both the bright-field and GFP settings (Figure 6).

*C. albicans treated with a fluorescent dye can autofluoresce, leading to a bright background (e.g., due to cell debris or the biofilm matrix) and bright cells dispersed over the viewing field. To take clear images, adjust the fluorescence excitation wavelength exposure time to reduce background. Dead cells will fluoresce green while those that do not fluoresce are live cells (Figure 6).*

16. Overlay bright-field and GFP images to observe the contrast between the live and dead cells. Take three representative images of each slide.

*Capture images with ~50 to 100 cells in the viewing field. Avoid including large clumps of cells in these images as the cells within them are often not distinguishable as separate entities and cannot be counted accurately.*





**Figure 6** A typical image generated by cell enumeration using fluorescence microscopy (Alternative Protocol 6) using fluorescein staining. The figure depicts both dead (**A**, green fluorescent) and live (**B**, dark non-fluorescent) *C. albicans* cells. Note also the faint background fluorescence which provides contrast to easily visualize the non-fluorescent live cells. Cells were imaged using standard GFP imaging settings on an EVOS-FL microscope. Scale bar = 100  $\mu\text{m}$ .

17. Using cell counting software, add the total number of live/dead cells in each captured image; there will be a total of nine images per data set.
18. Calculate mean, standard deviation, and statistical analyses (typically use a Student's unpaired two-tailed *t*-test assuming unequal variance) for each data set.

## BASIC PROTOCOL 9

### PHEROMONE-STIMULATED ASSAY

Recent work has shown that both mating type and pheromone signaling impact biofilm formation in *C. albicans*. *C. albicans* can exist as **a**,  $\alpha$ , or **a**/ $\alpha$  mating types, depending on the mating-type like (*MTL*) loci present. With the exception of this protocol, the biofilm assays described in this chapter pertain to biofilm formation of exclusively *MTL*-heterozygote **a**/ $\alpha$  *C. albicans* cells. Biofilm formation by **a**/ $\alpha$  cells has been reported to be more efficient than that by **a** or  $\alpha$  cells, indicating that genes present at the *MTL* loci affect biofilm development (Srikantha et al., 2012; Yi et al., 2011). In recent years, alternative *C. albicans* biofilms that are formed by *MTL*-homozygous white cells in response to pheromone released by cells of the opposite mating type have been characterized and these biofilms are often referred to as 'sexual biofilms' since they require the presence of mating pheromone (Daniels, Srikantha, Lockhart, Pujol, & Soll, 2006; Lin et al., 2013; Soll, 2014; Yi et al., 2011). The pheromone-stimulated assay presented below has been modified from the protocol initially reported by Lin et al. (2013). This assay is optimized for use with flat-bottomed, non-tissue culture-treated 12-well plates. It has been validated for use on clinically isolated white **a** and  $\alpha$  *C. albicans* strains, and can be compared to **a**/ $\alpha$  control strains used in the same assay. The pheromone-stimulated assay will allow researchers to test the biofilm forming abilities of *MTL*-homozygous strains in the presence of natural or synthetic pheromones of the opposite mating type or the biofilm abilities of co-cultures of both **a** and  $\alpha$  cells. In the 12-well format, three replicate wells are typically used per condition (blanks and controls included). Include blank wells containing medium alone, in addition to wells with cells not exposed to pheromone, as controls. We recommend that this assay be performed using  $\alpha$ -factor mating pheromone. Results using **a**-factor are much more variable due to its increased hydrophobicity.

## Materials

*C. albicans* strain(s) of interest

Spider medium, pH 7.2 (10 g/liter Difco brand nutrient broth, 10 g/liter mannitol, 4 g/liter K<sub>2</sub>HPO<sub>4</sub>, pH 7.2; store at room temperature)

Lee's medium, pH 6.8 (see recipe and Bedell & Sollet, 1979; Lin et al., 2013)

Dulbecco's phosphate-buffered saline (PBS; calcium and magnesium salt free), sterile filtered

Synthetic  $\alpha$  pheromone, sequence GFRLTNFGYFEPG, 13aa, 90% purity (LifeTein and see Bennett, Uhl, Miller, & Johnson, 2003)

YEPD agar plates (2% Bacto™ peptone, 2% dextrose, 1% yeast extract, 2% agar; store at room temperature)

12-well microtiter plates, non-tissue culture-treated (Falcon brand, manufacturer no. 351143)

Plate reader with optical density capabilities (600 nm) compatible with 12-well plate format (e.g., Tecan Infinite M1000 Pro or BioTek Epoch 2)

Roller drum for test tubes and/or shaker for flasks in a 25°C incubator

1. Streak *C. albicans* strains onto YEPD plates from glycerol stocks 2 to 5 days in advance of assay. Incubate plates at 30°C, 48hr.

*Do not use colonies that are >7 days old or store plates at 4°C as C. albicans can acquire aneuploidies under these conditions. New plates should be streaked from glycerol stocks rather than by re-streaking cells from existing plates.*

### **For small numbers of strains/large volumes**

- 2a. Inoculate a single colony in 25 ml Spider medium in a 250-ml flask (or 50 ml medium in a 500-ml flask) from a 2- to 5-day-old plate. Grow at 25°C with shaking overnight.

### **For large numbers of strains/small volumes**

- 2b. Inoculate a single colony in 4 ml Spider medium in a 20-ml test tube from a 2- to 5-day-old plate. Grow at 25°C with shaking or on a roller drum overnight.
3. After 16 to 18 hr, dilute overnight culture 1:20 or 1:40 and measure OD<sub>600</sub> using cuvettes and a spectrophotometer or 96-well plates and a plate reader.

*If assays will be set up >18 hr after overnight cultures are started, remove overnight cultures from the shaker or roller drum in the morning and allow to sit at room temperature until time of use (cultures should not sit for >4 hr at room temperature). Vortex cultures immediately before use. However, if possible, it is better to time the experiment such that the strains will be used prior to 18 hr of overnight growth.*

4. Following cell density determination for overnight cultures, add cells to a final OD<sub>600</sub> = 2.5 (or equivalent to  $\sim 5 \times 10^7$  cells/ml) in 1 ml Lee's medium.
5. Add synthetic *C. albicans*  $\alpha$  pheromone to each well at a final concentration of 10  $\mu$ g/ml (6.6  $\mu$ M). Distribute pheromone throughout the well by gently agitating the plate.
6. Incubate plate 24 hr at 25°C under static conditions.

*MTL-homozygous cells adhere very weakly before pheromone stimulation. As a result, this assay does not utilize the initial 90-min shaking adherence step and PBS wash employed by conventional biofilm assays. Pheromone-stimulated biofilms are generally less robust than conventional biofilms and thus are grown statically to minimize biofilm disruption.*

7. Gently aspirate medium from each well.

*Both the pheromone-stimulated and control biofilms will be very fragile. Try not to disturb the biofilm and leave a small layer of supernatant covering the biofilm.*

8. Gently wash biofilms with 1 ml PBS: Press pipet tip against the side of the well. Slowly add buffer while rotating the plate, ensuring even flow of PBS from all sides of the well. Then immediately aspirate PBS to prevent non-adhered cells from settling.

*Biofilms formed by MTL-homozygous *C. albicans* strains are generally very fragile. Be particularly careful during the wash step as entire biofilms can be washed away accidentally. Process one well at a time for this step.*

9. Measure OD<sub>600</sub> on a 12-well compatible plate reader.

*The number of reads per well may vary based on the instrument. We recommend taking the average density of five reads in each section of a five-by-five grid (i.e., twenty-five independent locations) in each well. Be sure that the wells are still moist for the reading as dry wells give inaccurate readings.*

10. Normalize data by subtracting the OD<sub>600</sub> reading of an average blank well (containing medium alone) from each experimental and control well. Normalize blank-subtracted OD<sub>600</sub> value of each experimental well to the mean blank-subtracted OD<sub>600</sub> for the relevant control wells.
11. Calculate the standard deviation and statistical analyses (typically use a Student's unpaired two-tailed *t*-test assuming unequal variance) for each normalized data set.

## **BASIC PROTOCOL 10**

### **TEMPORAL ASSESSMENT OF *C. ALBICANS* BIOFILM FORMATION USING A MICROFLUIDICS DEVICE**

This assay allows for the observation of biofilm formation as a biofilm develops over time using customizable conditions that mimic those of the host, such as those encountered in vascular catheters (Gulati, Ennis, Rodriguez, & Nobile, 2017; Lohse et al., 2017; Winter et al., 2016). This assay is highly customizable and alterations in temperature, media, and flow rates can be made to suit the needs of the researchers. The assay can be used to test the effects of antifungal compounds or to compare biofilm development between different strains. In our experience, there is good correlation between biofilm formation of mutant strains observed using this assay with that of the *in vivo* central venous rat catheter model (Andes et al., 2004; Gulati et al., 2017). Although this assay provides highly predictive data pertaining to biofilm development *in vivo*, it requires costly specialized equipment (the BioFlux 1000z) and has reduced throughput compared to many of the other methods described. We typically use two channels per strain or condition in a specialized 48-well plate (containing 24 channels), allowing for the evaluation of twelve conditions per plate (including blanks and controls).

#### **Materials**

RPMI-1640 medium with L-glutamine and 3-(*N*-morpholino)propanesulfonic acid (MOPS), without sodium bicarbonate, pH 7.0 (see recipe)

Compound(s) to be tested

70% isopropanol

Lens paper (VWR, cat. no. 5284-001)

Cuvettes and spectrophotometer or transparent 96-well plate and OD<sub>600</sub>-capable plate reader

BioFlux 1000z (Fluxion Biosciences)

48-well plate 0-20 dyne (Fluxion Biosciences, cat. no. 910-0047)

Montage software (Fluxion, Version 7.8.4.0)  
ImageJ software [National Institutes of Health (NIH)]

1. Initialize microfluidic device and set temperature to 37°C. Heat plate holder and plate to 37°C. Preheat RPMI-1640 medium at 37°C until the start of the experiment.

*Medium used for the experiment should be preheated to reduce the chances of air bubbles developing in the microfluidic device. Although the protocol listed here uses RPMI-1640 medium, other biofilm-inducing media can also be used.*

2. Clean interface plate and 0-20 dyne 48-well plate using 70% isopropanol and lens paper.

*The interface plate connects to the microfluidic system and the 48-well plate. Make sure to remove all liquids when cleaning to avoid blurry or dirty images and videos during the course of the experiment.*

3. Place plate on the microfluidic device and add 600 µl pre-heated RPMI-1640 to the desired inlet wells; we recommend at least two replicate wells per condition(s) to be tested.

*For a 12-hr experiment, 600 µl is sufficient. For longer experiments, add an additional 50 µl/hr, but do not exceed the maximum well volume of 1500 µl.*

4. Mount the clean interface onto the microfluidics plate, align it, and lock it in position using the levers to ensure that the system is airtight.

5. Using the computer interface for the program, flow medium from the inlet to the outlet wells at 1 dyne/cm<sup>2</sup>, 5 to 10 min to prime microfluidic channels.

*The microfluidic channels connect the inlet and outlet wells. Priming the channels with medium is essential to remove all air and avoid air bubbles that can interfere with the experiment.*

6. Following cell density determination for overnight cultures, dilute cells in a microcentrifuge tube to a final OD<sub>600</sub> = 0.5 (or equivalent to ~1 × 10<sup>7</sup> cells/ml) in 200 µl pre-warmed RPMI-1640 medium and preheat tube at 37°C until cells are added to the plate.

*C. albicans cells will begin forming hyphae as soon as they are resuspended in RPMI-1640 medium at 37°C. In order to capture images before filamentation begins, it is important to make the cell dilutions immediately prior to adding the cells to the microfluidic channels.*

7. Remove interface plate and add 50 µl of the cell cultures from step 6 to each outlet well.

*Take care to add the cells to the bottom of the well and not to introduce any air bubbles as these will disrupt the flow of medium and biofilm formation.*

8. Place interface plate back on the microfluidic plate and lock it in place.

9. Flow medium from the outlet to the inlet wells at 2 dyne/cm<sup>2</sup>, 3 sec to seed cells in the microfluidic channels.

*This step will flush C. albicans yeast cells into the viewing chamber, where the images will be acquired. It is important that the flow is stopped quickly to prevent contamination of the inlet wells supplying fresh medium. If these wells are contaminated, the channels supplying medium will be blocked with biofilms, confounding the results.*

10. Incubate cells without flow 20 min to allow cell adherence to the inner surface of the microfluidic channels.

11. During the 20-min adherence step, set up the stage positions on the computer software to allow for image acquisition during the experiment. Calibrate plate position and focus on each microfluidic channel to obtain the best image and save the settings.

*Three stage positions per channel are recommended, distributed evenly through the channel.*

12. Set acquisition for 145 total images at intervals of 5 min.

*We recommend capturing images at 50% Brightfield and 50% camera, with an exposure of 12 to 20 msec, 0.6 gain, and 20 MHz digitizer.*

*It is recommended that the experiment should be set up in dark conditions to avoid light interference during image capture.*

13. Begin image acquisition to capture the initial stages of biofilm development.
14. After the first round of images is captured, flush away weakly adhered cells by flowing medium from inlet to outlet wells at 1 dyne/cm<sup>2</sup>, 5 min.
15. Change flow rate to 0.5 dyne/cm<sup>2</sup> from inlet to outlet wells and leave experiment running 12 hr undisturbed and under dark conditions.
16. After 12 hr, use the computer software to view images and create a time-lapse video of biofilm formation (Figure 7 and Videos S1 and S2).


*See Gulati et al. (2017) for additional details on software acquisition settings.*

## **CONFOCAL IMAGING OF *C. ALBICANS* BIOFILMS**

Confocal scanning laser microscopy (CSLM) is a valuable tool in assaying biofilm formation as it allows for visualization of the biofilm architecture as well as measurement of the thickness of the biofilm (Nobile et al., 2012; Nobile & Mitchell, 2005). An array of different labeling methods are available to fluorescently stain biofilms. Polysaccharide-based staining methods—such as concanavalin A and Calcofluor-white, which bind to components of the fungal cell wall—are the most commonly used stains for *C. albicans* biofilms. Other stains that bind to DNA, proteins, and lipids are also available. The protocol described utilizes concanavalin A-Alexa Fluor 594 to stain *C. albicans* biofilms. Biofilms formed on microtiter plates or silicone squares can be visualized using this assay. For microtiter plates, we recommend that biofilms be formed on 6-well plates and that a blank well be included to monitor contamination during experimental setup. For CSLM imaging of both silicone squares and microtiter plates, two replicates per strain or condition are recommended. In this protocol, we describe imaging using an LSM 700 confocal microscope (Carl Zeiss). Other confocal microscopes may also be used, with the settings appropriate for each instrument. Many of the details in steps 9 through 13 will differ depending on the system used and as such these steps should be taken as general guidelines. This protocol is derived from Nobile et al. (2005) and Nobile et al. (2012).

### **Materials**

- C. albicans* strain(s) of interest
- RPMI-1640 medium with L-glutamine and 3-(*N*-morpholino)propanesulfonic acid (MOPS), without sodium bicarbonate, pH 7.0 (see recipe)
- Dulbecco's phosphate-buffered saline (PBS; calcium and magnesium salt free), sterile filtered
- Concanavalin A, Alexa Fluor 594 conjugate (Thermo Fisher Scientific, cat. no. C11253)



**Figure 7** Typical images for *C. albicans* wild-type (WT) and a known biofilm-defective mutant strain from the microfluidics assay (Basic Protocol 10). Using the BioFlux 1000Z, biofilms were grown for 12 hr post-adherence in Spider medium under dynamic flow ( $0.5 \text{ dyne/cm}^2$ ) at  $37^\circ\text{C}$ . Representative 0 (immediately post-adherence and initial wash), 4, 8, and 12 hr images (top to bottom) are shown for the WT strain (left column) and a known biofilm-defective mutant strain (right column). Scale bars =  $50 \mu\text{m}$ . The corresponding time-lapse videos of biofilm formation for these two strains were acquired at 15 frames/sec and are provided in Videos S1 and S2.

Transparent, sterile, flat-bottomed, non-tissue culture treated 6-well microtiter plates (Falcon brand, manufacturer no. 351146)  
 Serological pipets (VWR, cat. no. 89130-896)  
 ELMI DTS-4 shakers (or equivalent shaking incubators capable of holding 6-well plates, shaking at 200 rpm,  $37^\circ\text{C}$ )  
 Vacuum aspirator setup (to which 1000  $\mu\text{l}$  pipet tips can be attached)  
 Confocal microscope  
 ZEN software (Carl Zeiss)  
 ImageJ software [National Institutes of Health (NIH)]

1. Following cell density determination for overnight cultures, add cells to wells at a final  $\text{OD}_{600} = 0.5$  (equivalent to  $\sim 1 \times 10^7$  cells/ml) in 4 ml biofilm forming medium, such as RPMI-1640, Spider medium, or other medium of interest, in a 6-well plate.
2. Shake plate at  $37^\circ\text{C}$ , 90 min at 200 rpm in an ELMI incubator.

3. Aspirate medium and wash wells with 4 ml PBS.  
*Change pipet tips on aspirator between wells of different strains and/or conditions. Take care not to scrape the bottom of the well with the tip when aspirating.*
4. Aspirate PBS and add 4 ml fresh medium to each well.
5. Shake plate at 37°C, 24 hr at 200 rpm.
6. Add 20 µl concanavalin A (10 mg/ml; 50 µg/ml final concentration) directly to each well containing a *C. albicans* biofilm.  
*Take care not to disturb or disrupt the biofilm with the tip. After the addition of the dye, the plate should be kept in the dark throughout the experiment.*
7. Shake plate at 37°C, 60 min at 200 rpm in an ELMI incubator.
8. Place plate on the stage of the microscope inside a dark room. Take care not to disrupt the biofilm while moving the plate.
9. To visualize the concanavalin A, Alexa Fluor 594 conjugate (red fluorophore), use a 555-nm diode laser with a main beam splitter 405(T80/R20)/488/555/639 on the first channel and a pass band of 550 nm.
10. Visualize biofilms using a water-dipping 40× objective.  
*A water dipping 63× objective may also be used. It is essential that water-dipping objectives be used with this assay, as the oil-dipping objectives are not compatible for use with the biofilm medium.*
11. Obtain Z-Stacks at 652 × 652 pixels (160 × 160 µm) with a pinhole of 1 Airy unit, imaging every 0.5 µm; the stack thickness of a wild-type *C. albicans* biofilm should be between 240 to 290 µm.  
*We recommend obtaining 4 to 6 Z-stacks at different locations from each well. Make sure that Z-stacks begin from the bottom of the biofilm and not just at the limit of dye penetration.*
12. Analyze raw data, obtained as .czi files, using ZEN software.  
*The ZEN software is used to measure the depth/thickness of biofilms.*
13. In addition, analyze .czi files using the project stacks function of ImageJ to generate top-views and side-views of each stack (Figure 8).


**ALTERNATE  
PROTOCOL 7**

**CULTURING AND CONFOCAL IMAGING OF *C. ALBICANS*  
MIXED-SPECIES BIOFILMS**

Microbes in nature and in the human host are predominantly found as mixed-species biofilms, composed of many different microbes. This protocol allows for the assessment of the interactions of multiple species in the context of biofilms. Mixed-species biofilms formed with a standard optical density assay or silicone square assay can then be visualized using confocal microscopy. This allows for visualization of the architecture of the mixed species biofilm, determination of localization of each species within the biofilm, and measurement of the thickness or depth of the biofilm. Below we describe a protocol for assaying *C. albicans* and *E. coli* dual-species biofilms derived from Nobile & Mitchell, 2005 and Nobile et al., 2012, and Fox et al., 2014.

*Additional Materials (also see Basic Protocol 4 and Basic Protocol 11)*

SYTO 9 (5 mM in DMSO; Thermo Fisher Scientific, cat. no. S34854)



**Figure 8** Typical confocal scanning laser microscopy (CSLM) images of a *C. albicans* wild-type (WT) biofilm grown and imaged following Basic Protocol 11. Both the top-view (left) and side-view (right) of the WT biofilm are shown. Biofilms were grown in Spider medium in 6-well plates prior to staining with concanavalin A-Alexa Fluor 594 conjugate. A mixture of yeast and hyphal cells are visible in the top-view (left) but predominantly hyphae are visible in the side-view (right). The extracellular matrix is not visible in either image as it does not efficiently take up the stain. The stain does not penetrate to the base of the biofilm, hence the decrease in color deeper into the biofilm and the inability to visualize the yeast-form cells attached to the solid surface at the base of the substrate. Scale bars = 10  $\mu\text{m}$ .

1. Set up single and dual-species biofilms using a 6-well plate according to Basic Protocol 4, steps 1 through 9.
2. Add 20  $\mu\text{l}$  10 mg/ml (50  $\mu\text{g/ml}$  final concentration) of Concanavalin A-Alexa Fluor 594 directly to each well of *C. albicans* cultured alone or *C. albicans* cultured together with *E. coli* (or other bacterial species of interest).
3. Add 4  $\mu\text{l}$  Syto 9 nucleic acid stain (5  $\mu\text{M}$  final concentration) to wells of *E. coli* cultured alone or *C. albicans* cultured together with *E. coli*.

*Take care not to disturb or disrupt the biofilm with the tip. After the addition of the dye, the plate should be kept under dark conditions throughout the experiment.*

4. Shake plate at 37°C, 60 min at 200 rpm in an ELMI incubator.
5. Analyze biofilms using a confocal microscope as described in Basic Protocol 11, using a 488-nm diode laser on a second channel to visualize bacterial cells.

### HARVESTING CONDITIONED MEDIUM FROM *C. ALBICANS* BIOFILMS FOR PROTEOMIC ANALYSES

The biofilm conditioned medium assay offers a quick, relatively high-throughput way to harvest conditioned medium from biofilms for use with various proteomic analyses (Winter et al., 2016). This medium can then be processed and used in a variety of assays, e.g., multiplex substrate profiling by mass spectrometry (MSP-MS), and trypsin digests followed by shotgun proteomics (O'Donoghue et al., 2012).

We recommend the use of RPMI-1640 medium for this assay, however other media can also be used. We recommend avoiding media that contains serum and/or peptides (e.g., Spider medium, YEPD medium) as they prevent accurate quantification and normalization of samples and can interfere with subsequent analyses. We typically use eleven wells

**BASIC  
PROTOCOL 12**

**Gulati et al.**

**31 of 45**



(~44 ml, plus a 4-ml uninoculated control well) per strain, however a larger number of wells may be needed for samples that give low protein yields and/or for analyses that require more material. Depending on the comparisons to be made, it may be necessary to collect conditioned media from the same strain(s) grown under planktonic conditions. This protocol has been optimized for *C. albicans*, however it is possible to adapt it to work with other *Candida* species. When optimizing for other species, changes to the medium have the greatest effects on the assay.

### Materials

Dulbecco's phosphate-buffered saline (PBS; calcium and magnesium salt free), sterile filtered  
RPMI-1640 medium with L-glutamine and 3-(N-morpholino)propanesulfonic acid (MOPS), without sodium bicarbonate, pH 7.0 (see recipe)  
Liquid nitrogen  
Transparent, sterile, flat-bottomed, non-tissue culture treated 6-well microtiter plates (Falcon brand, manufacturer no. 351146)  
60-ml syringes (BD Falcon brand, manufacturer no. 309654)  
0.45- $\mu$ m PES membrane sterile syringe filters (Thermo Fisher Scientific, cat. no. 725-2545)  
Sterile 3-ml transfer pipets (BD Falcon brand, manufacturer no. 357575)  
50-ml screw cap tubes (e.g., Axygen Scientific SCT-50 ML-25-S or equivalent)  
ELMI DTS-4 shakers (or equivalent shaking incubators capable of holding 6-well plates, shaking at 200 rpm, 37°C)  
Vacuum aspirator setup (to which 200  $\mu$ l pipet tips can be attached)  
Electric pipettor (e.g., Drummond Pipet-Aid) and sterile 25 ml pipets  
Refrigerated centrifuge (capable of spinning 50 ml screw cap tubes)  
Cold room

1. Add 4 ml room temperature RPMI-1640 medium to each well of a 6-well plate.

*We typically use twelve wells (two plates) for the initial experiment with each strain for this assay.*

2. Following cell density determination for overnight cultures, add cells to wells at a final  $OD_{600} = 0.5$  (equivalent to  $\sim 1 \times 10^7$  cells/ml) in 4 ml RPMI-1640 medium. Do not inoculate first well on the last plate.

*We typically do not inoculate the first well on the last plate of a batch in order to control for contamination. When planktonic controls are needed, we inoculate to a final  $OD_{600} = 0.05$  in 25 ml in a 125-ml flask.*

3. Shake plate at 37°C, 90 min at 200 rpm in an ELMI incubator.

4. Aspirate medium. Try to remove any floating clumps of cells that are visible.

*Change pipet tips on aspirator between wells of different strains and/or plates. Take care not to scrape the bottom of the well with the tip when aspirating.*

5. Wash wells with 4 ml PBS. Aspirate PBS and try to remove any floating clumps of cells that are visible.

*Add PBS slowly to the side of the well, taking care to minimize disruption of cells on the surface of the well.*

6. Add 4 ml fresh medium to each well.

*Add medium slowly to the side of the well, taking care to minimize disruption of cells on the surface of the well.*

7. Shake plate at 37°C, 24 hr at 200 rpm in an ELMI incubator.
8. Pre-label 50-ml screw cap tubes. Chill screw cap tubes, 3-ml transfer pipets, syringes, and syringe filters on ice.

*We typically label and chill items used in the harvesting step 20 min before starting that step.*

9. Remove plates from ELMI shakers. Use sterile 3-ml plastic pipettors to scrape any cells off sides of the wells. Disrupt biofilm at the bottom of the well, collect both biofilm and liquid using the pipettor, and transfer both to the chilled, pre-labeled 50-ml tube. Keep 50-ml tube(s) on ice during this step.

*Change pipettors between plates. Materials from two plates should fit in a single 50-ml screw cap tube. We typically harvest from all wells of one plate before moving onto the next plate. Remove cells from the sides of all wells on a plate before collecting the liquid and biofilms from the wells. Depending on the strength of the biofilm, it may prove necessary to use the pipet to rigorously scrape the bottom of the well. There should not be visible cloudiness on the bottom of the well after the biofilms and conditioned medium are removed.*

10. Spin tubes 10 min at 3750 rpm in a chilled centrifuge. After spinning, put tubes on ice and transfer supernatant to new, pre-chilled tubes, taking care to minimize disruption of pellets. Keep new tubes with supernatants on ice.

*The robustness of the pellets can vary greatly between different strains, species, or conditions. Pour gently and do not worry about collecting all of the supernatant if this would risk transferring the cell pellet.*

11. Move to cold room (if not already there).
12. Pull plunger from 60-ml syringe and place it back in syringe packaging to avoid contamination. Attach 0.45- $\mu$ m PES filter to the syringe. Pour conditioned medium into syringe, reinsert plunger, and filter into fresh, pre-chilled 50-ml screw cap tube. Flash freeze tube in liquid nitrogen.

*Filter and flash freeze one tube before starting the next tube. Take care not to dislodge the filter from the syringe when filtering the medium.*

13. Store frozen samples at  $-80^{\circ}\text{C}$  pending further processing.

*Ideally processing and quantification of samples should be performed within a few days of harvesting. Details of further processing steps will vary depending on the requirements of subsequent proteomic analyses. When preparing conditioned medium for MSP-MS or trypsin digests followed by shotgun proteomics, we typically partially thaw sample in lukewarm tap water (until about  $\sim 50\%$  ice left) and then thaw the remainder on ice. We then use 10000 MWCO/15 ml spin units to concentrate samples down (one per tube) to  $\sim 0.5$  to 1 ml in a refrigerated centrifuge. Samples are then diluted to 15 ml with PBS and concentrated to  $\sim 1$  ml. Samples are then aliquoted, flash frozen with liquid nitrogen, and stored at  $-80^{\circ}\text{C}$ . Protein concentrations are determined using the Bradford assay in order to allow for normalization of input into later steps. If possible, use fresh samples for assays rather than ones that have been previously thawed and refrozen.*

## REAGENTS AND SOLUTIONS

*For common stock solutions, see Current Protocols, 2001.*

### **RPMI-1640 medium**

RPMI-1640 with L-glutamine, without sodium bicarbonate (MP Biomedicals, cat. no. 0910601) with 34.5 g/liter 3-(N-morpholino)propanesulfonic acid (MOPS). Adjust to pH 7.0 with sodium hydroxide.

Filter sterilize using 0.22- $\mu$ m filter. Store at 4°C per manufacturer's recommendations and protect from light.

### ***Lee's medium***

Final concentrations: 28.96% (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 1.15% MgSO<sub>4</sub>, 14.48% K<sub>2</sub>HPO<sub>4</sub>, 28.96% NaCl, 2.89% L-alanine, 7.53% L-leucine, 5.79% L-lysine, 0.57% L-methionine, 0.43% L-ornithine, 2.8% L-phenylalanine, 2.8% L-proline, 2.8% L-threonine, 0.49% L-arginine  
50 g/liter (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 2 g/liter MgSO<sub>4</sub>, 25 g/liter K<sub>2</sub>HPO<sub>4</sub>, 50 g/liter NaCl, 200 nM ZnSO<sub>4</sub>, 250 nM CuSO<sub>4</sub>, 1  $\mu$ M FeCl<sub>3</sub>, 1 mM MgCl<sub>2</sub>, biotin 10 mg/l and dextrose 1.25 g/liter. Adjust to pH 6.8. Store at room temperature for 1 to 2 months.

### ***Lethen broth***

Meat peptone 10 g/liter, beef extract 5 g/liter, polysorbate 80 5 g/liter, sodium chloride 5 g/liter, lecithin 0.7 g/liter. Store at room temperature for 1 to 2 months.

## **COMMENTARY**

### **Background Information**

#### ***Culturing***

We have presented common growth conditions that have worked across a number of strains and isolates in a large number of laboratory groups. YEPD medium is easy and inexpensive to make and can be stored for a long time. The growth conditions described consistently result in similar biofilms for a given strain.

#### ***Optical density assays***

These assays have been developed and optimized for medium to high throughput screening of deletion strain and potential antifungal compound libraries as a precursor to more time-consuming and/or cost-intensive assays (Lohse et al., 2017). The small volume, at least compared to many other biofilm assays, reduces the amount of test compound needed, which is particularly useful when testing expensive or limited compounds of interest. The dynamic ranges of these assays are also capable of detecting both increases and decreases in biofilm formation. They are also designed to minimize equipment requirements to one or more shaking incubators and a plate reader with optical density capabilities. Although we have presented several specific protocol variations (e.g., a 24-hr biofilm growth assay, a biofilm disruption assay, and several biofilm inhibition assays), the core protocol can be modified to study other stages of biofilm development (e.g., 4, 8, or 48 hr) or other metrics, such as synergistic drug reactions. Likewise, the related dispersal optical density assay can be modified to look at dispersal at other time points (Nobile et al., 2014).

The metric of biofilm formation varies between biofilm thickness (standard optical density assay), metabolic activity (XTT reduction assay), and cell viability (crystal violet assay) in the different assays; as such it may be informative to perform more than one variant. The standard optical density assay has the advantage of the minimal number of processing steps, both reducing the time involved and lessening the potential for inadvertent disruptions of the biofilm. It is important to remember, however, that there may be limitations to each of these approaches (e.g., not distinguishing between live and dead cells, missing metabolically inactive persister cells, and artifacts due to dyes failing to penetrate the full depth of the biofilm).

#### ***Mixed species assays***

*C. albicans* interacts with a wide range of bacterial species that occupy the same niches in humans (e.g., the oral cavity, skin, and the gastrointestinal and urogenital tracts; Gulati & Nobile, 2016; Nobile et al., 2015). *C. albicans* has been shown to form mixed species biofilms with several bacterial species both *in vitro* on abiotic surfaces and *in vivo* on biotic surfaces within the host (Harriott & Noverr, 2011; Morales et al., 2010). These interactions have significant clinical implications; for example, it has been observed that *C. albicans* biofilms can provide a protective anaerobic niche for pathogenic anaerobic bacteria, allowing these bacteria to grow in unexpected aerobic locations of the host when residing within biofilms with *C. albicans* (Fox et al., 2014). Likewise, *C. albicans* biofilms can protect bacteria from exposure to antibiotics (Harriott & Noverr, 2009, 2010). The modification of the standard optical

density assay to study *C. albicans* mixed-species biofilm formation with interacting bacterial partners offers opportunities to gain insights into cross-kingdom interactions occurring in the context of a biofilm. Furthermore, the base assay can be modified to study additional bacterial species as well as other aspects of biofilm development. Finally, these mixed species biofilm assays can be combined with the cell enumeration assays in order to evaluate the relative abundance and viability of each species within the biofilm.

#### **Silicone square assay**

*In vitro* *C. albicans* biofilm assays are most commonly performed on polystyrene microtiter plates. However, silicone is a common material used for implanted medical devices, such as catheters, heart pumps, and surgical reconstructive components. Silicone squares can be easily cut and customized for different experimental needs, and are thus quite versatile for non-high-throughput experimentation. The use of silicone squares is, therefore, a useful and potentially informative alternative substrate that can be worth testing in certain biofilm assays (Nobile et al., 2012; Nobile et al., 2005).

#### **Dry weight assay**

*In vitro* optical density assays offer a proxy for the number of cells and/or the amount of biomass present in a biofilm. The dry weight assay directly measures the amount of biomass that forms during the course of biofilm development and accounts for the presence of extracellular matrix in a way many assays cannot. Since this assay is low throughput, largely due to the vacuum filtration steps, it is used more often for validation of candidates initially identified by other methods.

#### **Cell enumeration assays**

Optical density assays offer, at best, crude proxies for the viability of cells in a biofilm, a limitation that can be of concern when screening potential antifungal compounds for microbicidal properties. An alternative assay, the cell enumeration assay, utilizes serial dilution and colony forming units to estimate the number of live cells present during biofilm formation. This approach has been utilized effectively in the cell adhesion assay to estimate the number of cells that adhere to seed a biofilm (Lohse et al., 2017; Winter et al., 2016). However, traditionally, it has been a challenge to obtain CFUs from a mature biofilm as biofilms are resilient to physical stresses, and as such

do not homogenize sufficiently, resulting in a substantial underestimation in the number of viable cells present. The addition of Lethen broth, which contains emulsifying components, allows the biofilm to break up and homogenize more effectively (without lysing the cells within the biofilm), resulting in a more accurate estimation of cell numbers. The alternate cell enumeration assay, which utilizes a fluorescent dye (e.g., fluorescein diacetate for *C. albicans*) to distinguish between live and dead cells, avoids the need to fully homogenize the biofilm by relying on the visual distinction between dead cells, which emit a bright green fluorescence and live cells, which are generally very dim (Figure 6). This version of the cell enumeration assay avoids the need for serial dilution, significantly reducing the workload and enabling the screening of more strains and/or conditions at one time. We note that this pair of assays is also helpful for looking at the relative abundance and viability of individual species within a mixed-species biofilm.

#### **Pheromone-stimulated assay**

This assay was originally established in Lin, Choi, and Bennett (2011) and optimized in Lin et al. (2013). It examines the ability of mating pheromone to stimulate biofilm formation of *C. albicans* *MTL*-homozygous strains, which would generally not result in robust biofilms in the absence of pheromone. The assay has been used to test *MTL*-homozygous strains using natural or synthetic pheromones of the opposite mating type, or co-cultures of both **a** cells and  **$\alpha$**  cells. In general, this assay works better using  **$\alpha$** -factor rather than **a**-factor due to  **$\alpha$** -factor's greater solubility in water. We note that even minor alterations to the steps of this assay can result in significant changes to the biofilms formed.

#### **Microfluidics assay**

Most of the *C. albicans* biofilm assays included in this article (e.g., the optical density, silicone square, and dry weight assays) assess biofilm formation at a single time point. To better assess *C. albicans* biofilms throughout their development, we recommend the real-time microfluidic assay using the BioFlux microfluidic device (Gulati et al., 2017). This assay is an effective tool to assess biofilm defects between different strains of *C. albicans* and to assess the effect of antifungal compounds on biofilm formation over time (Lohse et al., 2017; Winter et al., 2016). This assay is particularly useful to validate results from

high-throughput screens. In our experiences, the results from this assay have correlated well with *in vivo* experiments using the rat central venous catheter model (Andes et al., 2004; Winter et al., 2016).

### **Confocal imaging**

*C. albicans* forms highly structured biofilms that can be visualized by confocal scanning laser microscopy (CSLM). This technique allows for examination of the architecture of the biofilm as well as measurement of the biofilm thickness. CSLM has been utilized to visualize *C. albicans* biofilms in several studies (Nobile et al., 2012; Nobile et al., 2005), mostly utilizing concanavalin A conjugated with different Alexa Fluor™ dyes (most commonly Alexa Fluor 594). Recently, CSLM has also been used to study dual-species biofilm formation by staining fungal cells with concanavalin A dyes and bacterial strains with SYTO™ dyes (most commonly, SYTO™ 9; Fox et al., 2014). CSLM can be used to visualize biofilms formed on both microtiter plates and silicone squares.

### **Conditioned medium harvesting for proteomics**

The conditioned medium harvesting protocol presented here was developed to be as similar as possible to standard biofilm assays, while also satisfying the sample requirements of subsequent proteomic analyses. This protocol has worked effectively with a number of different strains and species (Winter et al., 2016).

## **Critical Parameters and Troubleshooting**

### **Culturing**

Standard microbiology sterile technique is critical for all protocols in this article using microbial cells and we commonly include one or more uninoculated negative controls to verify that medium has not been contaminated. Especially when working with a small number of strains, we recommend preparing multiple overnight cultures of each strain so that there are backup cultures in the event of individual culture contamination. When culturing *C. albicans*, it is critical to be aware that aneuploidies can be acquired if *C. albicans* is treated in ways commonly used for *Saccharomyces cerevisiae* and many bacterial species. In particular, it is important not to use colonies that are more than 7 days old, not to store plates at 4°C, and not to re-streak cells from existing plates. Likewise, strains should be used prior

to 18 hr of overnight growth. When using a spectrophotometer for the first time, it is important to first obtain CFU counts to correlate to optical density values. As overnight cultures will be very dense, it is important to dilute samples at least 20-fold to accurately determine the optical density. When working with mutant strains, slow growth phenotypes may result in reduced biofilm formation for reasons unrelated to the ability to form biofilms. As such, it is important to note any mutant(s) whose overnight cultures have significantly lower densities. In certain situations, these mutant strains could be allowed additional time to form a biofilm to fully assess biofilm development in a particular assay.

### **Optical density assays**

A wide range of experimental and even equipment-related factors can affect the output of the various optical density based biofilm assays. To help address this, we strongly recommend the inclusion of both positive and negative controls on every plate as well as controls for any solvents (e.g., DMSO, methanol) used to dissolve compounds. In general, we recommend validating the phenotypes of all mutants of interest using independently derived isolates and/or adback strains. Likewise, mutants should also be evaluated for other relevant phenotypes (e.g., adherence defects, filamentation defects) to classify possible reasons for defects in biofilm formation as well as to eliminate false positives (e.g., slow growth in general or in the medium being used). Attention should also be paid to the dynamic range of biofilm phenotypes observed when using a given medium; if the wild-type strain forms only a thin biofilm when using a certain medium, then it may be difficult to detect a reduction in biofilm formation with this medium. Alternatively, medium that supports robust biofilm formation may mask subtle defects of mutant strains. In these cases, altering the medium or specific nutrients (e.g., carbon sources) may be helpful.

There are a number of specific issues to be aware of when performing these assays. Variation in medium, the volume of medium, incubation times, and the way in which the medium is added or removed from wells can affect biofilm formation. As such, care should be taken to be as consistent as possible with each step. A related problem stems from the inadvertent physical disruption of biofilms when removing or adding liquid (e.g., the biofilm detaching during aspiration steps, during the

addition of liquid in the disruption optical density assay, or during replacement of the medium in the standard dispersal assay), which can range in scale from small disruptions at the edge of a well to the complete loss of a biofilm in a well. Any wells that are physically disrupted should be noted and excluded from subsequent analyses; the ability to exclude the occasional disrupted well is part of the reasoning behind having six to eight replicate wells of each condition. To avoid inadvertent disruption, it is important to remove or add medium slowly and to avoid touching the bottom of the well with pipet tips. Incomplete sealing of plates or peculiarities of specific incubators can result in position specific effects on biofilm formation; in general, we recommend repeating assays where this is noticeable. In order to reduce effects from prolonged exposure to air and the biofilm drying out, we recommend processing only one plate at a time, performing the disruption assay on only six to twelve wells at a time, and performing the final OD<sub>600</sub> measurements immediately after aspirating medium from the wells. When working with potential antifungal compounds, be sure to visually inspect the plate for the appearance of unusual colors and/or crystallization in wells as this may affect the OD<sub>600</sub> or OD<sub>492</sub> measurements. If such problems arise, it may be necessary to reduce the concentration of the compound being tested and/or try dissolving the compound in different solvents. Finally, it is important to remember that the output from the optical density, XTT reduction, or crystal violet assays may not scale linearly between different strain backgrounds, mutants, or species. As such, interesting phenotypes should be validated using additional methods, such as confocal microscopy or the microfluidics assay.

#### **Mixed species assays**

When setting up mixed-species biofilms, it is important to know the numbers of *C. albicans* and bacterial cells being seeded for the biofilm (for dual-species biofilms, we often begin with 1:1 ratios of each species). We recommend that an accurate correlation of OD<sub>600</sub> to cell density be established using a hemocytometer and/or the colony-forming units assay for each bacterial strain, as OD<sub>600</sub> can be influenced by cell shape and size. It is essential that the assay include single species biofilms in the same medium as the mixed-species biofilms, to aid in assessment of the mixed-species biofilms.

#### **Silicone square assay**

Including both positive and negative controls (e.g., wild-type and biofilm-defective mutant strains) as well as a blank control (uninoculated well) within each experiment is important to ensure proper experimental setup. Consistency in the treatment of samples is critical for reproducibility when performing this assay. It is necessary to note the top side of each silicone square, where the biofilm has been seeded; if the square gets flipped during manipulations, the biofilm will not form properly. Additional care must be taken when lifting the silicone squares during the wash step as well as when placing the silicone square into the fresh medium plate. For all steps involved, we recommend handling the silicone square only from the edges using pointed forceps to prevent disturbing the adhered cells. When placing the washed silicone squares into the fresh medium plate, be sure that the silicone squares are thoroughly submerged in the medium. If bubbles are present at the surface of the medium, push the silicone square below the bubbles to prevent the square from floating on top of the medium.

#### **Dry weight assay**

There are a number of steps which can result in variability when performing the dry weight assay. It is critical to recover all of the biofilm from the surface of the well when performing this assay. As such, second washes of the well may be necessary to collect biofilm fragments that escaped the first harvesting pass. Likewise, it is important to make sure that all of the biofilm is transferred from the pipet to the filtration membrane. Applying the biofilm to the center, rather than the edges, of the membrane is important as this reduces the chances of losing biofilm cells during later manipulations of the membrane. Finally, it is critical to dry all membranes thoroughly and uniformly in order to avoid variation due to residual moisture in the samples. Care should be taken to note any instances where dried biofilms have started to flake off from the membrane as this can affect dry weights if the biofilm fully detaches and is lost during weighing.

#### **Cell enumeration assays**

Colonies should only be counted once they have grown enough to be clearly visible to the naked eye and all plates should be counted at the same time if possible. As this assay is highly dependent on the homogenous disruption of biofilms and successful serial dilutions, there can be significant variability between

samples. We recommend at least four replicates per strain or condition and to plate multiple dilutions from each replicate. If problems persist with the serial dilution method, we suggest visualizing the homogenized sample on a microscope to ensure that there are no cell clumps remaining. If cell clumps remain despite vigorous homogenization efforts, it may be best to use the live/dead fluorescence version of the assay to circumvent this problem. As with the other assays described in this article, if the effect of an antifungal compound is being tested, it is important to include controls for the solvent (without the compound) to ensure that the measured effect is due to the compound and not the solvent concentration in the medium. It is also recommended to use low retention filter tips when making serial dilutions, as *C. albicans* biofilms tend to stick to the walls of pipet tips, which may result in inaccuracies during serial dilutions. For the visualization of live/dead cells, it is important to include both a positive control (e.g., a known antifungal drug) and a negative control (e.g., medium with no drug). The autofluorescent cells in the negative control can help in distinguishing live/dead cells in other samples as they provide a baseline of dead cell autofluorescence and the lack of fluorescence in live cells. When using this assay to distinguish between species in a mixed-species biofilm, it is important to make sure that only one species will grow on a given type of plate (e.g., *C. albicans* on plates containing antibiotics) or that it is possible to distinguish between colonies of the two species using a given medium. The final number of cells from the species used can vary by at least an order of magnitude, as such it may be necessary to increase or decrease the number of serial dilutions accordingly.

#### ***Pheromone-stimulated assay***

The reproducibility of this assay can be markedly improved by ensuring a uniform distribution of cells across the bottom of each well. To obtain this distribution, do not allow cells in the inoculum to settle. Gentle swirling of the plate immediately following addition of cells can aid in the even distribution of the cells across the surface of each well. The PBS wash step prior to measuring OD<sub>600</sub> is the most critical and sensitive step in this protocol as entire biofilms can be inadvertently washed away. To minimize this, carefully pipet PBS down the sides of the well, while slowly rotating the plate (the medium will tend to flow in evenly towards the center of the well). A layer of non-adhered cells may accumulate in the wells dur-

ing the wash step; a very gentle agitation after the PBS has been added can help disperse these non-adhered cells. Wells will rapidly dry out after medium has been aspirated, decreasing the effectiveness of the wash step. As such, we recommend aspirating medium and washing one well at a time in order to prevent wells from drying out.

#### ***Microfluidics assay***

Including controls within each batch of experiments, such as an untreated wild-type strain and controls for any solvents used, is important to ensure proper experimental setup. Setting the direction of flow between inlet and outlet wells is a critical component of performing this assay. Cells must flow from the outlet wells towards the inlet wells during seeding of the biofilm, and sterile medium must flow from the inlet wells towards the outlet wells during biofilm formation. This setup ensures that the sterile medium does not get contaminated and that the biofilm is formed within the microfluidic channel rather than inside the inlet wells. If the cells flow too fast or for a longer duration than indicated in the protocol, the cells could contaminate the inlet wells. If this occurs, the inlet channel may become blocked due to biofilm growth and biofilm development will not properly take place within the microfluidic channel. Other key components of performing this assay are to ensure that the camera focus and exposure time are accurately set for each channel and that the apparatus is not disturbed throughout the experiment. Additionally, this assay must take place in a dark room. Any alterations in lighting or vibrations can lead to blurry images and videos. Finally, prewarming the medium and cleaning the equipment are critical for preventing the development of air bubbles and for reducing issues during image acquisition.

#### ***Confocal imaging***

For CSLM, it is critical that the biofilms remain intact during the image acquisition process. The biofilms should be handled gently and kept in the dark until imaging is complete. Imaging should begin at the base of the biofilm and move towards the top. Care should be taken not to disturb the microscope or stage during the imaging process. We recommend at least two replicate wells per condition, and to obtain at least six Z-stack images per well from different areas of the well. When setting up Z-stack acquisition of images, care should be taken in obtaining the first image slice as the dye often cannot

penetrate deep into the biofilm due to its thickness. It is important to ensure that the base of the biofilm, rather than the base of dye penetration, has been reached before starting image acquisition. This can be achieved by visualizing the substrate (polystyrene plate or silicone square), which will appear quite distinct from the biofilm, and moving upwards until the substrate is out of view. Once the substrate is just out of view, this is a good position to set for the first image slice.

### **Conditioned medium harvesting for proteomics**

Although this assay is fairly robust in generating active protein mixtures, it is important to practice good sterile technique and to include at least one uninoculated well in each set of plates to act as a contamination control. As with all protein work, care should be taken to keep samples cold during the harvesting and processing steps to avoid the loss of protein activity due to enzymatic degradation and/or denaturation. Otherwise, the most important concern is to match growth times and medium conditions as much as possible between different strains and conditions. Changes to the medium may be necessary to improve biofilm formation for given strains; when making these changes, the effect on biofilm formation by all relevant strains should be noted. Likewise, care should be taken to avoid medium or reagents with large amounts of peptides or full length proteins, which will increase the background baseline (noise) of the assay. If the yield from the initial two plate assays is insufficient, it may be necessary to scale up to four or even eight plates for harvesting.

When harvesting, be sure to note whether control wells have signs of contamination. If contamination is observed, we recommend discarding that batch of conditioned medium and repeating the assay. Likewise, take care to note whether there are a noticeable number of free-floating planktonic cells in the well. The presence of at least some planktonic cells may be unavoidable, especially in the case of strains that form weaker biofilms. In these cases, the resulting profiling of the medium may represent a mix of the biofilm and planktonic cell types.

## **Anticipating Results**

### **Culturing**

With some exceptions, colonies of a strain should have consistent phenotypes when streaked onto a plate. Cultures should look vis-

ibly saturated after overnight growth and the uninoculated negative control cultures should remain clear. In most cases, there should not be visible flocculants (clumps of cells) in the overnight cultures. A normal cell density measured by OD<sub>600</sub> for an overnight culture is around 10 to 12, although this can vary for certain mutant strains and based on the spectrophotometer (and path length) used to obtain the measurement.

### **Optical density assays**

Plates should be examined visually for noticeable phenotypes and abnormal results in addition to the more quantitative OD<sub>600</sub> and OD<sub>492</sub> measurements from a plate reader. There should be a visible biofilm on the bottom of positive control wells while the negative control wells should remain clear (Figure 2). Contaminants in the negative control wells may appear as discreet spots on an otherwise clear well or as a cloudy haze throughout the well. Ideally, the density of biofilms should appear relatively uniform across a given well and between replicate wells of a given condition. Mutants or compounds that reduce biofilm formation should result in less dense or even clear wells. Wells where the addition of the compound of interest leaves a noticeable color, visible precipitation, or crystallization should be noted as this may affect OD<sub>600</sub> or OD<sub>492</sub> measurements. Position specific effects (e.g., noticeably lower density for all wells on one side of a plate) should be noted. Expected OD<sub>600</sub> values will vary depending on the strain, medium, duration of growth, and even the plate reader; wild-type densities in the 0.2 to 0.6 range are not unexpected.

### **Mixed species assays**

The results from this assay are conceptually similar to those from the optical density assay and its variants. It is possible that different bacterial species may not form robust single species biofilms in the medium of choice (e.g., BHI-FBS) under aerobic conditions for a variety of reasons (e.g., nutritional or oxygen requirements). Despite this, the measurements from these single species control biofilms are important for contrasting against the dual-species biofilms formed by the same species in combination with *C. albicans* under the identical conditions. As the OD<sub>600</sub> results do not distinguish cell density measurements between the two species, it will often be necessary to pair this assay with the cell enumeration assays.



### ***Silicone square assay***

Silicone squares should be examined visually to ensure expected biofilm formation and expected turbidity of the surrounding medium. Positive control strains capable of forming normal biofilms should have clear surrounding medium with the silicone substrate completely covered by the biofilm (Figure 4). Negative control squares unable to form a biofilm should have dense turbidity of the surrounding medium with the silicone substrate devoid of cells. Blank (uninoculated) wells should be completely clear both in the medium as well as on the silicone square (Figure 4).

### ***Dry weight assay***

There should be a visible collection of biofilm cells on the dried membranes from the positive control wells and no visible cells on the dried membranes from the negative control wells. Although major biofilm defects will be clearly detectable using this assay, more minor biofilm defects will be less apparent or sometimes undetectable. In instances where a minor biofilm defect is anticipated, other biofilm assays should be used to assess the biofilm phenotype.

### ***Cell enumeration assays***

When counting colonies, the plates examined within 24 hr should contain single colonies (not lawns) and should be easily distinguished from one another. If no colonies are observed after 24 hr, plates should be examined daily for 48 to 72 hr, which will allow for the scoring of slow growing mutant strains. When using this assay with mixed-species biofilms, the results should be visually distinct colonies for each species or colonies of single species that only grow on specific selective medium for that species.

### ***Pheromone-stimulated assay***

The results of this assay are conceptually similar to the optical density assays. Wells that have not been treated with pheromone should appear mostly clear after the PBS wash, with a thin layer of yeast cells remaining on the bottom of the well. Wells treated with pheromone will form a fragile biofilm of varying thickness depending on strain background. Strains that respond strongly to pheromone will form a biofilm that coats the entire well; very few, if any, non-adherent cells should be washed away in those cases. Strains with a moderate response to pheromone will form a biofilm with upper layers that detach around the edges. Weak responders to pheromone will

have an entire layer of non-adherent cells that comes off the well (oftentimes as a sheet) during the wash, with only the most adherent cells remaining. Strains that do not respond to pheromone will be washed off the well in small clumps of cells rather than sheets of biofilm.

### ***Microfluidics assay***

At the completion of the experiment, the microfluidic plate can be examined visually. The inlet wells should have clear medium and the outlet wells should have an obvious biofilm (in the case of a biofilm-forming strain). Additionally, most of the medium should be in the outlet wells, with less than 100  $\mu$ l present in the inlet wells. If any of the inlet wells still have noticeable amounts of medium, it is possible that the channel became blocked during the experiment due to an air bubble. If the inlet wells have a biofilm, however, this indicates that the cells reached and contaminated the inlet well during the biofilm seeding step; the corresponding microfluidic channel should not be included in further analysis. The images generated at the end of the experiment can be analyzed using the Montage software. Images should not be blurry and the viewing field should not be obstructed by debris (Figure 7). During the generation of a time-lapse video using the Montage software, ensure that there were no air bubbles evident during the experiment. The video should begin by showing adhered yeast-form cells, followed by the formation of hyphae within 30 min (in the case of wild-type strains in biofilm-inducing medium), followed by the formation of a thick biofilm covering most of the field of view at the end of the 12-hr experiment (Video S1).

### ***Confocal imaging***

Wild-type *C. albicans* biofilms are highly structured, containing intercalating hyphae visible throughout the Z-stack of images (Figure 8). Usually yeast-form cells can be seen interspersed between the hyphae (Figure 8). Depending on the medium and substrate used, the thickness of a typical wild-type biofilm after 24 hr can range from  $\sim$ 100 to 400  $\mu$ m (and can be much thicker if grown for 48 to 60 hr). Mutant strains or strains exposed to antifungal agents may have a much thinner biofilm. When imaging mutant strains, it is important to note if particular cell types are evident in excess or in significantly reduced amounts relative to the wild-type strain.

### ***Conditioned medium harvesting for proteomics***

It is not possible to visually evaluate the harvested conditioned medium at the end of the harvesting or processing steps. The exact concentration of the processed medium will depend on the volume of medium harvested and the final volume after the processing step; we typically get a yield of around 500  $\mu\text{l}$  at 100 to 400  $\mu\text{g/ml}$  from 50 ml of conditioned medium. This material is consistently active in *in vitro* protease cleavage assays and at least several dozen distinct proteins are anticipated to be detected when shotgun proteomics are performed on a trypsin digest of the processed medium.

### **Time Considerations**

#### ***Culturing***

The time required for streaking plates from frozen ( $-80^{\circ}\text{C}$ ) stocks, starting overnight cultures, and determining the densities of cultures will vary depending on the number of strains being used. Small numbers ( $<10$ ) of strains will take only a few minutes for each of these steps.

#### ***Optical density assays***

The most time-consuming steps of these assays are often the preparation of solutions (when screening compounds for anti-biofilm activity) and or the dilutions of strains (when screening different isolates or mutant strains). The length of these steps will depend greatly on the number of dilutions being made. Transferring these solutions to deep well 96- or 384-well plates can take between 10 min (for 96-well) and 30 min (for 384-well) for each plate. Once solutions and cells are in a suitable format or in cases where medium and/or cells can be added from a sterile basin, starting a plate takes  $\sim 15$  min. Washing and adding fresh medium to the plate takes between 15 min (for 96-well) and 30 min (for 384-well) for each plate. Other than the dilution step, the overall setup process should be expected to take roughly 2.5 hr with active work required for less than half that time. Once familiar with this technique, it is possible to stagger the setup of two plates by offsetting them  $\sim 45$  min and to set up multiple staggered sets of plates in a day. The disruption step, when applicable, will take between 30 to 45 min per plate. Aspirating the medium and reading the density of wells will commonly take between 15 min (for 96-well) and 30 min (for 384-well) for each plate. Solutions for the XTT reduction assay will take  $\sim 15$  min to prepare. Aspirat-

ing the existing medium and adding the XTT solution to wells will take 15 min (for 96-well) to 30 min (for 384-well) plates. After the incubation step, which may vary in length (assume 30 min for initial assays), reading the plate will take  $\sim 5$  min. The crystal violet assay takes around 3 to 4 hr, although much of this time is taken up by the 45 min drying, staining, and destaining steps. Harvesting medium for the dispersal assays takes  $\sim 10$  min per plate. These times are approximate and are based on familiarity with the assay; steps may take considerably longer when performing the assay for the first time.

#### ***Mixed species assays***

This assay will take longer at first for any new bacterial species due to the need to correlate the  $\text{OD}_{600}$  with actual cell numbers using a hemocytometer, which can take  $\sim 30$  min per strain. Alternatively, the  $\text{OD}_{600}$ -cell number correlation can be determined in advance using the colony-forming units assay. Once the  $\text{OD}_{600}$  to cell number correlation is established, this assay should take roughly the same amount of time as the standard optical density assay. Depending on the specific bacterial strain(s), however, there may be additional preparatory work required for culturing the specific species.

#### ***Silicone square assay***

Determining the cell density, placing the silicone squares into the wells, and adding medium and cells to wells takes  $\sim 20$  min for one 12-well plate. Fresh 12-well plates with PBS and RPMI-1640 can be prepared during the 90-min adhesion step. The PBS wash step takes up the most time and care, about 1 to 3 min per silicone square or 12 to 36 min per 12-well plate.

#### ***Dry weight assay***

It will take  $\sim 5$  to 10 min to set up and  $\sim 10$  min to wash each 6-well plate. Depending on the nature and robustness of the biofilms, it will take between 5 to 10 min to disrupt and homogenize each plate. It will take several minutes to set up the filter device (Figure 5), filter the biofilm, and remove the membrane for each sample. As such, processing times of at least 15 min can be expected for each 6-well plate.

#### ***Cell enumeration assays***

The serial dilution/plate counting version of this assay takes significantly longer than the fluorescent dye based version. The most time-consuming steps of the former version of

the assay are the serial dilutions and the plating of multiple dilutions for all replicates. A dilution series for four replicates of a given strain will generally take ~30 to 40 min. Once colonies have grown on the plates, it can take 1 to 5 min to count each plate depending on colony density. For the fluorescent dye-based protocol, the imaging step is the most time-intensive step of this assay. It takes ~15 min to acquire all representative images for each slide. Unlike the serial diluting/plate counting version, however, cell counting software, such as imageJ (NIH), make the analysis of each slide relatively fast (5 to 10 min per image).

#### ***Pheromone-stimulated assay***

Most of the time requirements for this assay are similar to those for the related standard optical density assay. Seeding biofilms will take ~1 hr, which is longer than the standard optical density assay, due to the pheromone addition step. The final washing step can take ~45 min due to the care required while washing. Analysis with a spectrophotometer takes ~10 min for each plate.

#### ***Microfluidics assay***

This assay can be set up in <1 hr once the researcher is familiar with all the steps and the software. The experiment takes 12 hr to complete and the analysis, including the generation of time-lapse videos, can take 2 to 4 hr depending on the number of strains and conditions used in the experiment.

#### ***Confocal imaging***

This assay is time intensive, with the bulk of time spent on image acquisition using the confocal microscope. It takes ~15 to 20 min to obtain one set of Z-stack images, so each well can take ~3 hr to image.

#### ***Conditioned medium harvesting for proteomics***

Determining cell density, adding medium, and adding cells will take ~15 to 20 min for a set of four 6-well plates (24 wells total). The wash step takes ~5 to 10 min per plate, or ~30 min for a set of four 6-well plates. As such, the full setup process takes ~2.5 hr with active work for less than half of the time. Physically harvesting medium from wells takes ~5 to 10 min per plate, or ~20 to 45 min for a set of four 6-well plates depending in part on the effort needed to disrupt the biofilm. Coupled with the spin and filtering steps, the full harvesting process will take between 30 to 60 min for a set of four 6-well plates. These steps be-

come faster with practice; once familiar with the assay, it is possible to stagger the setup, washing, and processing of two sets of four 6-well plates. Subsequent processing steps will vary in length; the method we use takes ~2 hr for thawing and at least 3 hr for the concentration steps.

#### **Acknowledgements**

We thank all members of the Nobile lab and BioSynthesis team for helpful discussions on biofilm protocols and optimization procedures. This work was supported by National Institutes of Health (NIH) grants R21AI125801 (to C.J.N.), R35GM124594 (to C.J.N.), and R43AI131710 (to M.B.L.). C.L.E. was partially supported by a fellowship from the Sonoma County Mycological Association. D.L.R. was supported by a doctoral fellowship from the University of California Institute for Mexico and the United States (UC-MEXUS) and Consejo Nacional de Ciencia y Tecnologia (CONACYT).

#### **Literature Cited**

- Andes, D., Nett, J., Oschel, P., Albrecht, R., Marchillo, K., & Pitula, A. (2004). Development and characterization of an in vivo central venous catheter *Candida albicans* biofilm model. *Infection and Immunity*, 72(10), 6023–6031. doi: 10.1128/IAI.72.10.6023-6031.2004.
- Baillie, G. S., & Douglas, L. J. (1999). Role of dimorphism in the development of *Candida albicans* biofilms. *Journal of Medical Microbiology*, 48(7), 671–679. doi: 10.1099/0022-2615-48-7-671.
- Bedell, G. W., & Soll, D. R. (1979). Effects of low concentrations of zinc on the growth and dimorphism of *Candida albicans*: Evidence for zinc-resistant and -sensitive pathways for mycelium formation. *Infection and Immunity*, 26(1), 348–354.
- Bennett, R. J., Uhl, M. A., Miller, M. G., & Johnson, A. D. (2003). Identification and characterization of a *Candida albicans* mating pheromone. *Molecular and Cellular Biology*, 23(22), 8189–8201. doi: 10.1128/MCB.23.22.8189-8201.2003.
- Calderone, R. A., & Fonzi, W. A. (2001). Virulence factors of *Candida albicans*. *Trends in Microbiology*, 9(7), 327–335. doi: 10.1016/S0966-842X(01)02094-7.
- Chandra, J., Kuhn, D. M., Mukherjee, P. K., Hoyer, L. L., McCormick, T., & Ghannoum, M. A. (2001). Biofilm formation by the fungal pathogen *Candida albicans*: Development, architecture, and drug resistance. *Journal of Bacteriology*, 183(18), 5385–5394. doi: 10.1128/JB.183.18.5385-5394.2001.
- Current Protocols. (2001). Commonly used reagents. *Current Protocols in*

- Microbiology*, A.2A.1–A.2A.15. doi: 10.1002/9780471729259.mca02as00.
- Daniels, K. J., Srikantha, T., Lockhart, S. R., Pujol, C., & Soll, D. R. (2006). Opaque cells signal white cells to form biofilms in *Candida albicans*. *EMBO Journal*, 25, 2240–2252. doi: 10.1038/sj.emboj.7601099.
- Douglas, L. J. (2002). Medical importance of biofilms in *Candida* infections. *Revista Iberoamericana De Micologia*, 19(3), 139–143.
- Douglas, L. J. (2003). *Candida* biofilms and their role in infection. *Trends in Microbiology*, 11(1), 30–36. doi: 10.1016/S0966-842X(02)00002-1.
- Finkel, J. S., Xu, W., Huang, D., Hill, E. M., Desai, J. V., Woolford, C. A., . . . Mitchell, A. P. (2012). Portrait of *Candida albicans* adherence regulators. *PLoS Pathogens*, 8(2), e1002525. doi: 10.1371/journal.ppat.1002525.
- Fox, E. P., Bui, C. K., Nett, J. E., Hartooni, N., Mui, M. C., Andes, D. R., . . . Johnson, A. D. (2015). An expanded regulatory network temporally controls *Candida albicans* biofilm formation. *Molecular Microbiology*, 96(6), 1226–1239. doi: 10.1111/mmi.13002.
- Fox, E. P., Cowley, E. S., Nobile, C. J., Hartooni, N., Newman, D. K., & Johnson, A. D. (2014). Anaerobic bacteria grow within *Candida albicans* biofilms and induce biofilm formation in suspension cultures. *Current Biology*, 24(20), 2411–2416. doi: 10.1016/j.cub.2014.08.057.
- Gulati, M., Ennis, C. L., Rodriguez, D. L., & Nobile, C. J. (2017). Visualization of biofilm formation in *Candida albicans* using an automated microfluidic device. *Journal of Visualized Experiments*, (130). doi: 10.3791/56743.
- Gulati, M., & Nobile, C. J. (2016). *Candida albicans* biofilms: Development, regulation, and molecular mechanisms. *Microbes and Infection*, 18(5), 310–321. doi: 10.1016/j.micinf.2016.01.002.
- Harriott, M. M., & Noverr, M. C. (2009). *Candida albicans* and *Staphylococcus aureus* form polymicrobial biofilms: Effects on antimicrobial resistance. *Antimicrobial Agents and Chemotherapy*, 53(9), 3914–3922. doi: 10.1128/AAC.00657-09.
- Harriott, M. M., & Noverr, M. C. (2010). Ability of *Candida albicans* mutants to induce *Staphylococcus aureus* vancomycin resistance during polymicrobial biofilm formation. *Antimicrobial Agents and Chemotherapy*, 54(9), 3746–3755. doi: 10.1128/AAC.00573-10.
- Harriott, M. M., & Noverr, M. C. (2011). Importance of *Candida*-bacterial polymicrobial biofilms in disease. *Trends in Microbiology*, 19(11), 557–563. doi: 10.1016/j.tim.2011.07.004.
- Hawser, S. P., & Douglas, L. J. (1994). Biofilm formation by *Candida* species on the surface of catheter materials in vitro. *Infection and Immunity*, 62(3), 915–921.
- Jin, Y., Yip, H. K., Samaranyake, Y. H., Yau, J. Y., & Samaranyake, L. P. (2003). Biofilm-forming ability of *Candida albicans* is unlikely to contribute to high levels of oral yeast carriage in cases of human immunodeficiency virus infection. *Journal of Clinical Microbiology*, 41(7), 2961–2967. doi: 10.1128/JCM.41.7.2961-2967.2003.
- Kim, J., & Sudbery, P. (2011). *Candida albicans*, a major human fungal pathogen. *Journal of Microbiology*, 49(2), 171–177. doi: 10.1007/s12275-011-1064-7.
- Krom, B. P., Cohen, J. B., McElhaney-Feser, G., Busscher, H. J., van der Mei, H. C., & Cihlar, R. L. (2009). Conditions for optimal *Candida* biofilm development in microtiter plates. *Methods in Molecular Biology*, 499, 55–62. doi: 10.1007/978-1-60327-151-6\_7.
- Krom, B. P., Cohen, J. B., McElhaney Feser, G. E., & Cihlar, R. L. (2007). Optimized candidal biofilm microtiter assay. *Journal of Microbiological Methods*, 68(2), 421–423. doi: 10.1016/j.mimet.2006.08.003.
- Krom, B. P., & Willems, H. M. (2016). In vitro models for *Candida* biofilm development. *Methods in Molecular Biology*, 1356, 95–105. doi: 10.1007/978-1-4939-3052-4\_8.
- Kuhn, D. M., Balkis, M., Chandra, J., Mukherjee, P. K., & Ghannoum, M. A. (2003). Uses and limitations of the XTT assay in studies of *Candida* growth and metabolism. *Journal of Clinical Microbiology*, 41(1), 506–508. doi: 10.1128/JCM.41.1.506-508.2003.
- Kullberg, B. J., & Oude Lashof, A. M. (2002). Epidemiology of opportunistic invasive mycoses. *European Journal of Medical Research*, 7(5), 183–191.
- Kumamoto, C. A. (2002). *Candida* biofilms. *Current Opinion in Microbiology*, 5(6), 608–611. doi: 10.1016/S1369-5274(02)00371-5.
- LaFleur, M. D., Kumamoto, C. A., & Lewis, K. (2006). *Candida albicans* biofilms produce antifungal-tolerant persister cells. *Antimicrobial Agents and Chemotherapy*, 50(11), 3839–3846. doi: 10.1128/AAC.00684-06.
- LaFleur, M. D., Lucumi, E., Napper, A. D., Diamond, S. L., & Lewis, K. (2011). Novel high-throughput screen against *Candida albicans* identifies antifungal potentiators and agents effective against biofilms. *Journal of Antimicrobial Chemotherapy*, 66(4), 820–826. doi: 10.1093/jac/dkq530.
- Lafleur, M. D., Sun, L., Lister, I., Keating, J., Nantel, A., Long, L., . . . Lewis, K. (2013). Potentiation of azole antifungals by 2-adamantanamine. *Antimicrobial Agents and Chemotherapy*, 57(8), 3585–3592. doi: 10.1128/AAC.00294-13.
- Lin, C. H., Choi, A., & Bennett, R. J. (2011). Defining pheromone-receptor signaling in *Candida albicans* and related asexual *Candida* species. *Molecular Biology of the Cell*, 22(24), 4918–4930. doi: 10.1091/mbc.E11-09-0749.
- Lin, C. H., Kabrawala, S., Fox, E. P., Nobile, C. J., Johnson, A. D., & Bennett, R. J. (2013). Genetic control of conventional and pheromone-stimulated biofilm

- formation in *Candida albicans*. *PLoS Pathogens*, 9(4), e1003305. doi: 10.1371/journal.ppat.1003305.
- Lohse, M. B., Gulati, M., Johnson, A. D., & Nobile, C. J. (2018). Development and regulation of single- and multi-species *Candida albicans* biofilms. *Nature Reviews Microbiology*, 16(1), 19–31. doi: 10.1038/nrmicro.2017.107.
- Lohse, M. B., Gulati, M., Valle Arevalo, A., Fishburn, A., Johnson, A. D., & Nobile, C. J. (2017). Assessment and optimizations of *Candida albicans* in vitro biofilm assays. *Antimicrobial Agents and Chemotherapy*, 61(5). doi: 10.1128/AAC.02749-16.
- Lopez-Ribot, J. L. (2005). *Candida albicans* biofilms: More than filamentation. *Current Biology*, 15(12), R453–455. doi: 10.1016/j.cub.2005.06.020.
- Morales, D. K., & Hogan, D. A. (2010). *Candida albicans* interactions with bacteria in the context of human health and disease. *PLoS Pathogens*, 6(4), e1000886. doi: 10.1371/journal.ppat.1000886.
- Nailis, H., Kucharikova, S., Ricicova, M., Van Dijck, P., Deforce, D., Nelis, H., & Coenye, T. (2010). Real-time PCR expression profiling of genes encoding potential virulence factors in *Candida albicans* biofilms: Identification of model-dependent and -independent gene expression. *BMC Microbiology*, 10, 114. doi: 10.1186/1471-2180-10-114.
- Nett, J. E., Cain, M. T., Crawford, K., & Andes, D. R. (2011). Optimizing a *Candida* biofilm microtiter plate model for measurement of antifungal susceptibility by tetrazolium salt assay. *Journal of Clinical Microbiology*, 49(4), 1426–1433. doi: 10.1128/JCM.02273-10.
- Nobile, C. J., Andes, D. R., Nett, J. E., Smith, F. J., Yue, F., Phan, Q. T., ... Mitchell, A. P. (2006). Critical role of Bcr1-dependent adhesins in *C. albicans* biofilm formation in vitro and in vivo. *PLoS Pathogens*, 2(7), e63. doi: 10.1371/journal.ppat.0020063.
- Nobile, C. J., Fox, E. P., Hartooni, N., Mitchell, K. F., Hnisz, D., Andes, D. R., ... Johnson, A. D. (2014). A histone deacetylase complex mediates biofilm dispersal and drug resistance in *Candida albicans*. *MBio*, 5(3), e01201–01214. doi: 10.1128/mBio.01201-14.
- Nobile, C. J., Fox, E. P., Nett, J. E., Sorrells, T. R., Mitrovich, Q. M., Hernday, A. D., ... Johnson, A. D. (2012). A recently evolved transcriptional network controls biofilm development in *Candida albicans*. *Cell*, 148(1-2), 126–138. doi: 10.1016/j.cell.2011.10.048.
- Nobile, C. J., & Johnson, A. D. (2015). *Candida albicans* biofilms and human disease. *Annual Reviews Microbiology*, 69, 71–92. doi: 10.1146/annurev-micro-091014-104330.
- Nobile, C. J., & Mitchell, A. P. (2005). Regulation of cell-surface genes and biofilm formation by the *C. albicans* transcription factor Bcr1p. *Current Biology*, 15(12), 1150–1155. doi: 10.1016/j.cub.2005.05.047.
- Nobile, C. J., Nett, J. E., Andes, D. R., & Mitchell, A. P. (2006). Function of *Candida albicans* adhesin Hwp1 in biofilm formation. *Eukaryotic Cell*, 5(10), 1604–1610. doi: 10.1128/EC.00194-06.
- Nobile, C. J., Nett, J. E., Hernday, A. D., Homann, O. R., Deneault, J. S., Nantel, A., ... Mitchell, A. P. (2009). Biofilm matrix regulation by *Candida albicans* Zap1. *PLoS Biology*, 7(6), e1000133. doi: 10.1371/journal.pbio.1000133.
- Norice, C. T., Smith, F. J., Jr., Solis, N., Filler, S. G., & Mitchell, A. P. (2007). Requirement for *Candida albicans* Sun41 in biofilm formation and virulence. *Eukaryotic Cell*, 6(11), 2046–2055. doi: 10.1128/EC.00314-07.
- O'Donoghue, A. J., Eroy-Reveles, A. A., Knudsen, G. M., Ingram, J., Zhou, M., Statnekov, J. B., ... Craik, C. S. (2012). Global identification of peptidase specificity by multiplex substrate profiling. *Nature Methods*, 9(11), 1095–1100. doi: 10.1038/nmeth.2182.
- Pierce, C. G., Chaturvedi, A. K., Lazzell, A. L., Powell, A. T., Saville, S. P., McHardy, S. F., & Lopez-Ribot, J. L. (2015). A novel small molecule inhibitor of *Candida albicans* biofilm formation, filamentation and virulence with low potential for the development of resistance. *NPJ Biofilms Microbiomes*, 1, pii: 15012. doi: 10.1038/npjbiofilms.2015.12.
- Pierce, C. G., Saville, S. P., & Lopez-Ribot, J. L. (2014). High-content phenotypic screenings to identify inhibitors of *Candida albicans* biofilm formation and filamentation. *Pathogens and Disease*, 70(3), 423–431. doi: 10.1111/2049-632X.12161.
- Ramage, G., Vande Walle, K., Wickes, B. L., & Lopez-Ribot, J. L. (2001). Standardized method for in vitro antifungal susceptibility testing of *Candida albicans* biofilms. *Antimicrobial Agents and Chemotherapy*, 45(9), 2475–2479. doi: 10.1128/AAC.45.9.2475-2479.2001.
- Richard, M. L., Nobile, C. J., Bruno, V. M., & Mitchell, A. P. (2005). *Candida albicans* biofilm-defective mutants. *Eukaryotic Cell*, 4(8), 1493–1502. doi: 10.1128/EC.4.8.1493-1502.2005.
- Soll, D. R. (2014). The evolution of alternative biofilms in an opportunistic fungal pathogen: An explanation for how new signal transduction pathways may evolve. *Infection, Genetics and Evolution*, 22, 235–243. doi: 10.1016/j.meegid.2013.07.013.
- Srikantha, T., Daniels, K. J., Pujol, C., Sahni, N., Yi, S., & Soll, D. R. (2012). Nonsex genes in the mating type locus of *Candida albicans* play roles in  $\alpha$  biofilm formation, including impermeability and fluconazole resistance. *PLoS Pathogens*, 8(1), e1002476. doi: 10.1371/journal.ppat.1002476.
- Tournu, H., & Van Dijck, P. (2012). *Candida* biofilms and the host: Models and new concepts for eradication. *International Journal of Microbiology*, 2012, 845352. doi: 10.1155/2012/845352.

- Uppuluri, P., Chaturvedi, A. K., Srinivasan, A., Banerjee, M., Ramasubramaniam, A. K., Kohler, J. R., ... Lopez-Ribot, J. L. (2010). Dispersion as an important step in the *Candida albicans* biofilm developmental cycle. *PLoS Pathogens*, *6*(3), e1000828. doi: 10.1371/journal.ppat.1000828.
- Uppuluri, P., Pierce, C. G., Thomas, D. P., Bubeck, S. S., Saville, S. P., & Lopez-Ribot, J. L. (2010). The transcriptional regulator Nrg1p controls *Candida albicans* biofilm formation and dispersion. *Eukaryotic Cell*, *9*(10), 1531–1537. doi: 10.1128/EC.00111-10.
- Winter, M. B., Salcedo, E. C., Lohse, M. B., Hartooni, N., Gulati, M., Sanchez, H., ... Nobile, C. J. (2016). Global Identification of Biofilm-Specific Proteolysis in *Candida albicans*. *MBio*, *7*(5). doi: 10.1128/mBio.01514-16.
- Yi, S., Sahni, N., Daniels, K. J., Lu, K. L., Srikantha, T., Huang, G., ... Soll, D. R. (2011). Alternative mating type configurations (a/alpha versus a/a or alpha/alpha) of *Candida albicans* result in alternative biofilms regulated by different pathways. *PLoS Biology*, *9*(8), e1001117. doi: 10.1371/journal.pbio.1001117.