

Draft Genome Sequence of *Leucobacter* sp. Strain UCD-THU (Phylum *Actinobacteria*)

Hannah E. Holland-Moritz, Dakota R. Bevans, Jenna M. Lang, Aaron E. Darling,* Jonathan A. Eisen, David A. Coil

University of California Davis Genome Center, Davis, California, USA

* Present address: Aaron E. Darling, University of Technology Sydney, Ultimo, NSW, Australia.

H.E.H.-M. and D.R.B. contributed equally to this work.

Here we present the draft genome of *Leucobacter* sp. strain UCD-THU. The genome contains 3,317,267 bp in 11 scaffolds. This strain was isolated from a residential toilet as part of an undergraduate project to sequence reference genomes of microbes from the built environment.

Received 21 April 2013 Accepted 26 April 2013 Published 6 June 2013

Citation Holland-Moritz HE, Bevans DR, Lang JM, Darling AE, Eisen JA, Coil DA. 2013. Draft genome sequence of *Leucobacter* sp. strain UCD-THU (phylum *Actinobacteria*). *Genome Announc.* 1(3):e00325-13. doi:10.1128/genomeA.00325-13.

Copyright © 2013 Holland-Moritz et al. This is an open-access article distributed under the terms of the [Creative Commons Attribution 3.0 Unported license](https://creativecommons.org/licenses/by/3.0/).

Address correspondence to Jonathan A. Eisen, jaeisen@ucdavis.edu.

Members of the *Leucobacter* genus are yellow-pigmented, Gram-positive, aerobic, nonmotile, non-spore-forming, irregular rod-shaped bacteria (1). They can be found in a variety of environments, including nematodes (2), waste treatment facilities (3), sediments receiving chromium-contaminated water (4), and cow dung (5).

Leucobacter sp. strain UCD-THU was isolated from a residential toilet in Davis, California, in an effort to provide reference genomes of microbes that are found in the built environment (6, 7). Biofilm scrapings were incubated overnight in Luria broth (LB) at 37°C, then plated out on LB agar. Single colonies were picked for serial dilution streaking and the organism was identified by Sanger sequencing of the 16S rRNA gene after PCR amplification using the 1391R and 27F primers. Genomic DNA was extracted using a Wizard genomic DNA purification kit (Promega) from fresh overnight cultures. Illumina paired-end libraries were then made from sonicated DNA using a TruSeq DNA sample prep v2 kit (Illumina). Fragments between 300 and 600 bp were selected using a Pippin Prep DNA size selection system (Sage Science).

A total of 3,121,654 paired-end reads were generated on an Illumina MiSeq, at a read length of 160 bp. Quality trimming and error correction of the reads resulted in 2,942,884 high-quality reads. All sequence processing and assembly were performed using the a5 assembly pipeline. This pipeline automates the processes of data cleaning, error correction, contig assembly, scaffolding, and quality control. The assembly produced 103 contigs, contained in 11 scaffolds (minimum, 5,587 bp; maximum, 699,082 bp; N_{50} , 443,492 bp). During scaffolding, some contigs were merged based on short overlaps and read-pair information, yielding a final collection of 30 contigs in 11 scaffolds that were submitted to GenBank. This resulted in a final assembly of 3,317,267 bp with a GC content of 70% and an overall coverage estimate of 280×. Completeness of the genome was assessed using the PhyloSift software (A. E. Darling, G. Jospin, E. Lowe, E. Matsen, H. Bik, and J. A. Eisen, submitted), which searches for a

list of 40 highly conserved, single-copy marker genes (D. Wu, G. Jospin, and J. A. Eisen, in preparation), of which all 40 were found in this assembly.

Automated annotation was performed using the RAST annotation server (8). *Leucobacter* sp. strain UCD-THU contains 3,018 predicted coding sequences and 52 predicted RNAs. A full-length (1,488-bp) 16S sequence was obtained from this annotation and was used to attempt to identify the species of *Leucobacter*. The top 250 BLAST hits were aligned with MUSCLE (9), which was used to construct a phylogenetic tree with FastTree 2 (10) (doi:10.6084/m9.figshare.674598). *Leucobacter* sp. strain UCD-THU is most closely related to a single isolate, *L. chironomi*, but is more than 99% identical to several other *Leucobacter* species; therefore, we are unable to assign a species name to this isolate.

The genome sequence of only one other *Leucobacter* species has been published—that of *L. chromiirestis* (11); this species also possesses a 16S rRNA gene sequence that is 99% identical to *Leucobacter* sp. strain UCD-THU.

Nucleotide sequence accession numbers. This whole-genome shotgun project has been deposited at DDBJ/EMBL/GenBank under the accession number [APJM000000000](https://www.ncbi.nlm.nih.gov/nuccore/APJM000000000). The version described in this paper is the first version, [APJM010000000](https://www.ncbi.nlm.nih.gov/nuccore/APJM010000000). Illumina reads are available at doi:10.6084/m9.figshare.201772.

ACKNOWLEDGMENTS

Illumina sequencing was performed at the DNA Technologies Core facility in the Genome Center at the UC Davis, Davis, California.

This work was funded by a grant to J.A.E. from the Alfred P. Sloan Foundation as part of their “Microbiology of the Built Environment” program.

REFERENCES

1. Takeuchi M, Weiss N, Schumann P, Yokota A. 1996. *Leucobacter komagatae* gen. nov., sp. nov., a new aerobic gram-positive, nonsporulating rod with 2,4-diaminobutyric acid in the cell wall. *Int. J. Syst. Bacteriol.* 46:967–971.

2. Somvanshi VS, Lang E, Schumann P, Pukall R, Kroppenstedt RM, Ganguly S, Stackebrandt E. 2007. *Leucobacter iarius* sp. nov., in the family *Microbacteriaceae*. *Int. J. Syst. Evol. Microbiol.* 57:682–686.
3. Halpern M, Shakéd T, Pukall R, Schumann P. 2009. *Leucobacter chironomi* sp. nov., a chromate-resistant bacterium isolated from a chironomid egg mass. *Int. J. Syst. Evol. Microbiol.* 59:665–670.
4. Morais PV, Paulo C, Francisco R, Branco R, Paula Chung A, da Costa MS. 2006. *Leucobacter luti* sp. nov., and *Leucobacter alluvii* sp. nov., two new species of the genus *Leucobacter* isolated under chromium stress. *Syst. Appl. Microbiol.* 29:414–421.
5. Weon HY, Anandham R, Tamura T, Hamada M, Kim SJ, Kim YS, Suzuki K, Kwon SW. 2012. *Leucobacter denitrificans* sp. nov., isolated from cow dung. *J. Microbiol.* 50:161–165.
6. Bendiks ZA, Lang JM, Darling AE, Eisen JA, Coil DA. 2013. Draft genome sequence of *Microbacterium* sp. strain UCD-TDU (phylum *Actinobacteria*). *Genome Announc.* 1(2):e0012013-13. doi:10.1128/genomeA.00120-13.
7. Lo JR, Lang JM, Darling AE, Eisen JA, Coil DA. 2013. Draft genome sequence of an actinobacterium, *Brachybacterium muris* strain UCD-AY4. *Genome Announc.* 1(2):e00086-13. doi:10.1128/genomeA.00086-13.
8. Aziz RK, Bartels D, Best AA, DeJongh M, Disz T, Edwards RA, Formsma K, Gerdes S, Glass EM, Kubal M, Meyer F, Olsen GJ, Olson R, Osterman AL, Overbeek RA, McNeil LK, Paarmann D, Paczian T, Parrello B, Pusch GD, Reich C, Stevens R, Vassieva O, Vonstein V, Wilke A, Zagnitko O. 2008. The RAST server: rapid annotations using subsystems technology. *BMC Genomics* 9:75.
9. Edgar RC. 2004. MUSCLE: a multiple sequence alignment method with reduced time and space complexity. *BMC Bioinformatics* 5:113.
10. Price MN, Dehal PS, Arkin AP. 2010. FastTree AAP 2—approximately maximum-likelihood trees for large alignments. *PLoS One* 5:e9490.
11. Sturm G, Buchta K, Kurz T, Rensing SA, Gescher J. 2012. Draft genome sequence of *Leucobacter chromiirensistens*, an extremely chromium-tolerant strain. *J. Bacteriol.* 194:540–541.