UC Office of the President

Recent Work

Title

Multiscale imaging of HIV-1 transmission in humanized mice

Permalink

https://escholarship.org/uc/item/273484zb

Authors

Kieffer, Collin Ladinsky, Mark S Bjorkman, Pamela J

Publication Date

2017-09-07

DOI

10.1089/AID.2017.0142

Peer reviewed

This paper has been peer-reviewed and accepted for publication, but has yet to undergo copyediting and proof correction. The final published version may differ from this proof

Multiscale imaging of HIV-1 transmission in humanized mice
(DOI: 10.1089/AID.2017.0142)

© Mary Ann Liebert, Inc.

DOI: 10.1089/AID.2017.0142

Multiscale imaging of HIV-1 transmission in humanized mice

Kieffer, C.¹, Ladinsky, M.S.¹, Bjorkman, P.J.¹

 Department of Biology and Biological Engineering, California Institute of Technology, Pasadena, CA.

Running Title: Multiscale Imaging of HIV-1 in Humanized Mice.

KEYWORDS: fluorescence microscopy, animal models, HIV/SIV Pathogenesis, electron microscopy, fluorescence microscopy, HIV transmission, Humanized mice, Tissue Clearing, Imaging

Corresponding author:

P.J. Bjorkman

California Institute of Technology

MC 114-96, Department of Biology and Bioengineering

1200 E. California Blvd.

Pasadena, CA 91125 USA

626-395-8350

bjorkman@caltech.edu


paper has been peer-reviewed and accepted for publication, but has yet to undergo copyediting and proof correction. The final published version may differ from this proof

Multiscale imaging of HIV-1 transmission in humanized mice
 (DOI: 10.1089/AID.2017.0142)

References:

- 1) Kieffer, C, Ladinsky, MS, Ninh, A, Galimidi, RP, Bjorkman, PJ. Longitudinal imaging of HIV-1 spread in humanized mice with parallel 3D immunofluorescence and electron tomography. Elife 2017; 6 pii: e23282.
- 2) Yang, B, Treweek, JB, Kulkarni, RP, Deverman, BE, Chen, CK, Lubeck, E, Shah, S, Cai, L, Gradinaru, V. Single-Cell Phenotyping within Transparent Intact Tissue through Whole-Body Clearing. Cell 2014; 158(4):945-58.
- 3) Treweek, J.B.; Chan, K.Y.; Flytzanis, N.C.; Yang, B.; Deverman, B.E.; Greenbaum, A.; Lignell, A.; Xiao, C.; Cai, L.; Ladinsky, M.S.; Bjorkman, P.J.; Fowlkes, C.C.; Gradinaru, V. Whole-body tissue stabilization and selective extractions via tissue-hydrogel hybrids for high-resolution intact circuit mapping and phenotyping. *Nat Prot*, 2015;10(11):1860-96.
- 4) Ladinsky, MS, Kieffer, C, Olson, G, Deruaz, M, Vrbanac, V, <u>Tager</u>, AM, <u>Kwon</u>, DS, <u>Bjorkman</u>, PJ. Electron tomography of HIV-1 infection in gut-associated lymphoid tissue. *PLoS Pathogens* 2014;10(1):e1003899.

paper has been peer-reviewed and accepted for publication, but has yet to undergo copyediting and proof correction. The final published version may differ from this proof



Despite ongoing efforts to characterize HIV-1 pathogenesis, there remain numerous gaps in our understanding of the biological mechanisms of HIV-1 dissemination in tissues. To visualize the structural details of HIV-1 transmission we developed a multiscale imaging platform to longitudinally observe virus dissemination at the resolution of individual infected cells and virions in lymphoid tissues from mice with humanized immune systems (hu-mice). A) Representative bioluminescence image showing the distribution of human Tcells within a hu-mouse (image courtesy of Rachel Galimidi). Immunodeficient mice (NOD.Cg-Prkd^{scid}II2rq^{tm1WjI}/SzJ; NOD/SCID/IL2Ry-/-; NSG) were injected with stimulated human peripheral blood mononuclear cells to rapidly generate a hu-mouse predominantly reconstituted with activated T-cells [1]. Hu-mice were infected with R5-tropic HIV-1 (YU-2), sacrificed at specific times post-infection, and lymphoid tissues including gut-associated lymphoid tissue (GALT), spleen, and reproductive tract were harvested and fixed for analysis by light and electron microscopy (EM). B) Fixed GALT samples (top) were cleared using the passive CLARITY technique [2, 3] to render tissues visually transparent (bottom) and entire samples were immunostained to identify human cells and viral proteins [1]. The removal of opaque biomolecules during clearing enhances light penetration into samples and refractive index matching of tissues, imaging media, and objectives minimizes optical distortions during imaging. These methods were applied to visualize large volumes (mm³cm³) of intact tissues with single cell resolution using confocal [1] and light sheet fluorescence microscopy (LSFM). Dashed box indicates approximate size of region imaged in subsequent panel. C) 3D-immunoflluorescence (3D-IF) image of GALT from an HIV-1infected hu-mouse captured by LSFM, which utilizes cylindrical lenses to generate thin sheets of light to optically section whole tissues. This method enabled the acquisition of

4

AIDS Research and Human Retroviruses
Multiscale imaging of HIV-1 transmission in humanized mice (doi: 10.1089/AID.2017.0142)
This article has been peer-reviewed and accepted for publication, but has yet to undergo copyediting and proof correction. The final published version may differ from this proof.

Downloaded by Rensselaer Polytechnic Institute NERL from online.liebertpub.com at 09/11/17. For personal use only.

AIDS Research and Human Retroviruses

Multiscale imaging of HIV-1 transmission in humanized mice
 (DOI: 10.1089/AID.2017.0142)

has been peer-reviewed and accepted for publication, but has yet to undergo copyediting and proof correction. The final published version may differ from this

paper

proof.

3D-IF images with single cell resolution and delineated the spatial relationships of HIVinfected cells and target cells within complex tissue environments. Human CD3+ T-cells (magenta) were distributed throughout GALT, and HIV-1 p24 staining (green) preferentially localized to colon crypts at the base of villi (Dashed box; also indicates approximate size of region imaged in ensuing panel). D) Adjacent regions of HIV-1-infected GALT from hu-mice were prepared for EM by a hybrid method involving high-pressure freezing and freezesubstitution fixation [4, 1]. Thin-section EM overview of a comparable region of colon crypts identified an area (orange box) that contained HIV-1 virions. E) Electron tomography (ET) of the virus containing region from D) revealed 3D ultrastructural information about the localization of individual virions and provided sufficient resolution to distinguish both mature and immature free-virions (white arrows) and nascent budding virions (black arrow). Note the orders of magnitude in change of scale between the various imaging methods used to analyze the same sample. We are currently applying these techniques to image additional tissues of interest (i.e. lymph nodes, thymus, brain, and bone marrow) from highly relevant animal models of HIV-infection and human clinical samples. Taken together, these efforts highlight the utility of multiscale imaging for understanding both the global pathology and structural details of HIV-1 transmission within lymphoid tissues, and provide new insights for combating HIV spread throughout the body.